

Product datasheet for **SC329527**

TGF beta induced factor 2 (TGIF2) (NM_001199513) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: TGF beta induced factor 2 (TGIF2) (NM_001199513) Human Untagged Clone
Tag: Tag Free
Symbol: TGIF2
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC329527 representing NM_001199513.
Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
ATGTCGGACAGTGATCTAGGTGAGGACGAAGGCCTCCTCTCCCTGGCGGGCAAAGGAAGCGCAGGGGG  
AACCTGCCAAGGAGTCGGTGAAGATCCTCCGGGACTGGCTGTACTTGACCCGCTACAACGCCTACCCC  
TCAGAGCAGGAGAAGCTGAGCCTTTCTGGACAGACCAACCTGTCAGTGCTGCAAATATGTAAGTGGTTC  
ATCAATGCCCGGGCGGGCTTCTCCAGACATGCTTCGGAAGGATGGCAAAGACCCTAATCAGTTTACC  
ATTTCCCGCCGCGGGGTAAAGGCCTCAGATGTGGCCCTCCCCGTGGCAGCAGCCCCTCAGTGCTGGCT  
GTGTCTGTCCAGCCCCCAATGTGCTCCTCCTGTCTGTGTGCTCCATGCCGCTTCACTCAGGCCAG  
GGGAAAAGCCAGCAGCCCTTTCCACGTGGGAGCTGGAGTCTCCAAGCCCTGGTGACCCCTGGT  
AGCACACTTACTCTGCTGACCAGGGCTGAGGCTGGAAGCCACAGGTGGACTCTTCAACACGCCACCA  
CCCACACCCAGAGCAGGACAAAGAGGACTTCAGCAGCTTCCAGCTGCTGGTGGAGGTGGCGCTACAG  
AGGGCTGCTGAGATGGAGCTTCAAGAGCAGGACCCATCACTCCCATTACTGCACACTCCCATCCCT  
TTAGTCTCTGAAAATCCCAGTAG
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:



ACCN: NM_001199513

Insert Size: 714 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001199513.1</u>
RefSeq Size:	3328 bp
RefSeq ORF:	714 bp
Locus ID:	60436
UniProt ID:	<u>Q9GZN2</u>
Cytogenetics:	20q11.23
Protein Families:	Transcription Factors
MW:	25.9 kDa
Gene Summary:	<p>The protein encoded by this gene is a DNA-binding homeobox protein and a transcriptional repressor, which appears to repress transcription by recruiting histone deacetylases to TGF beta-responsive genes. This gene is amplified and over-expressed in some ovarian cancers. Alternative splicing results in multiple transcript variants. A related pseudogene has been identified on chromosome 1. Read-through transcription also exists between this gene and the neighboring downstream C20orf24 (chromosome 20 open reading frame 24) gene. [provided by RefSeq, Dec 2010]</p> <p>Transcript Variant: This variant (3) differs in the 5' UTR compared to variant 1. All variants (1-4) encode the same protein.</p>