

## Product datasheet for **SC329573**

### DAP3 (NM\_001199849) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** DAP3 (NM\_001199849) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** DAP3  
**Synonyms:** bMRP-10; DAP-3; MRP-S29; MRPS29; S29mt  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC329573 representing NM\_001199849.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

ATGATGCTGAAAGGAATAACAAGGCTTATCTCTAGGATCCATAAGTTGGACCCTGGGCGTTTTTACAC
ATGGGGACCCAGGCTCGCCAAAGCATTGCTGCTCACCTAGATAACCAGGTTCCAGTTGAGAGTCCGAGA
GCTATTTCCCGCACCAATGAGAATGACCCGGCCAAGCATGGGGATCAGCACGAGGGTCAGCACTACAAC
ATCTCCCCCAGGATTTGGAGACTGTATTTCCCATGGCCTTCTCCTCGCTTTGTGATGCAGGTGAAG
ACATTCAGTGAAGCTTGCTGATGGTAAGGAAACCAGCCCTAGAACTTCTGCATTACCTGAAAAACACC
AGTTTTGCTTATCCAGCTATACGATATCTTCTGTATGGAGAGAAGGGAACAGGAAAAACCTAAGTCTT
TGCCATGTTATTCATTTCTGTGCAAAACAGGACTGGCTGATACTACATATTCAGATGCTCATCTTTGG
GTGAAAAATTGTCGGGATCTTCTGCAGTCCAGCTACAACAACAGCGCTTTGATCAACCTTTAGAGGCT
TCAACCTGGCTGAAGAATTTCAAACACTACAAATGAGCGCTTCTGAACCAGATAAAAAGTTCAAGAGAAG
TATGTCTGGAATAAGAGAGAAAGCACTGAGAAAGGGAGTCTCTGGGAGAAGTGGTTGAACAGGGCATA
ACACGGGTGAGGAACGCCACAGATGCAGTTGGAATTGTGCTGAAAGAGCTAAAGAGGCAAAGTCTTTG
GGTATGTTTACCTCCTAGTGGCCGTGGATGGAATCAATGCTCTTTGGGAAGAACCCTCTGAAAAAGA
GAAGATAAAAGCCCGATTGCCCCGAGGAATTAGCACTTGTTCACAACTGAGGAAAATGATGAAAAAT
GATTGGCATGGAGCGCCATTGTGTCGGCTTTGAGCCAGACTGGGCTCTCTTTAAGCCCGGAAAGCC
TATCTGCCCCAGGAGTTGCTGGGAAAGGAAGATTTGATGCCCTGGATCCCTTTATCCCATCCTGGTT
TCCAACATAAACCAGGAATTTGAAAGTTGTATTAGTATTATTTGAAAAACAATTGGCTTCAACAT
GAGAAAGCTCCTACAGAAGAAGGAAAAAAGAGCTGCTGTTCTAAGTAACGCGAACCCCTCGTGCTG
GAGCGGCACTGTGCCTACCTCTAA
  
```

**Restriction Sites:** SgfI-MluI



[View online »](#)

**Plasmid Map:**


**ACCN:** NM\_001199849

**Insert Size:** 1197 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001199849.1](#)

**RefSeq Size:** 2161 bp

**RefSeq ORF:** 1197 bp

<b>Locus ID:</b>	7818
<b>UniProt ID:</b>	<a href="#">P51398</a>
<b>Cytogenetics:</b>	1q22
<b>Protein Families:</b>	Druggable Genome
<b>MW:</b>	45.6 kDa
<b>Gene Summary:</b>	<p>Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 28S subunit protein that also participates in apoptotic pathways which are initiated by tumor necrosis factor-<math>\alpha</math>, Fas ligand, and gamma interferon. This protein potentially binds ATP/GTP and might be a functional partner of the mitoribosomal protein S27. Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. Pseudogenes corresponding to this gene are found on chromosomes 1q and 2q. [provided by RefSeq, Dec 2010]</p> <p>Transcript Variant: This variant (3) has an additional exon in the 5' UTR, as compared to variant 2. Variants 1, 2 and 3 encode the same isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>