

Product datasheet for **SC329605**

NDUFS1 (NM_001199981) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NDUFS1 (NM_001199981) Human Untagged Clone
Tag:	Tag Free
Symbol:	NDUFS1
Synonyms:	CI-75k; CI-75Kd; MC1DN5; PRO1304
Vector:	pCMV6-Entry (PS100001)




[View online »](#)

Fully Sequenced ORF: >SC329605 representing NM_001199981.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGTTAAGGATACCTGTAAGAAAGGCCTTAGTAGGCCTTCTAAGTCTCCTAAAGGATGTGTTCGAACA
ACTGCCACAGCAGCAAGCAACTTGATTGAAGTATTTGTTGATGGTCAGTCTGTCATGGTGGAAACCGGA
ACGACCGTCCTCCAAGTTGTAGCTGCTTGTGCCATGCCAGTAATGAAGGGTTGGAATATCCTAACAAC
TCAGAAAAATCCAAAAAGCCAGGGAAGGTGTGATGGAGTTCTTATTAGCAAATCACCCATTGGACTGT
CCTATTTGTGACCAGGGAGGTGAATGTGATCTGCAGGACCAGTCCATGATGTTTGGAAATGATAGGAGC
CGATTTTTAGAGGGGAAGCGTGTGCTGTGGAAGACAAGAACATTGGGCCATTGGTAAAGACCATCATGACA
AGATGTATACAGTGTACTCGCTGCATCAGGTTTCAAGTGAGATTGCAGGAGTAGATGATTTGGGAACA
ACAGGCAGAGGAAATGATATGCAAGTTGGCACATACATTGAAAAGATGTTTCATGTCTGAAGTGTCTGGG
AATATCATTGATATCTGCCCTGTAGGTGCCCTAACCTCTAAGCCCTATGCCTTTACTGCCCGCCTTGG
GAAACAAGAAAGACAGAATCCATTGATGTAATGGATGCGGTTGGAAGTAATATTGTGTTAGCACAAGA
ACTGGAGAAGTGATGAGGATTTGCCACGTATGCATGAGGACATCAATGAAGAGTGGATCTCTGATAAA
ACCAGATTTGCCATGATGGGCTAAAACGTCAAAGACTTACCGAGCCAATGGTCAGAAATGAAAAAGGG
CTTTTAACCTATACTTCTTGGGAGGATGCGCTCTCTCGCGTAGCTGGAATGTTGCAGAGTTTTCAAGGC
AAAGATGTGGCAGCAATTGCAGGTGGCTTGGTGGATGCTGAAGCCCTGGTAGCTCTCAAAGATTTGCTT
AATAGAGTGGACTCTGACACCTTATGCACTGAAGAGGTCTTCCCACTGCAGGAGCTGGCACAGATTTG
CGTTCCAATTAATCTTAACTACAATTGCTGGTGTGGAAGAGGCAGATGTTGTTCTTCTGGTTGGT
ACAAACCACGTTTTGAGGCACCACTGTTAATGCTAGAATTCGAAAGAGCTGGCTGCATAATGACTTA
AAAGTGGCCCTTATAGGCAGTCCAGTGGACCTCACTTACACATATGACCACCTGGGAGACTCCCCAAA
ATTCTTCAAGACATTGCTTCGGGAAGCCATCCATTTAGCCAGTCTTAAAGGAAGCTAAAAACCAATG
GTGGTTTTAGGCAGTTCTGCACTCAAAGAAATGATGGAGCAGCAATTCTTGCAGCTGTTTCTAGCATT
GCACAAAAGATTCGGATGACTAGTGGTGTACTGGTATTGGAAGTTATGAATATCCTTCATAGGATT
GCAAGTCAAGTAGCTGCTTTGGACCTTGCTATAAGCCTGGGGTGAAGCAATTCGGAAGAACCTCCC
AAGGTGCTGTTTCTCCTGGGAGCAGATGGAGGTTGTATCACACGACAGGATTTGCCAAAGGATTGTTTC
ATTATTTATCAAGGACATCATGGTGTGTTGGGGCTCCCATAGCTGATGTTATTCTCCAGGAGCTGCT
TACACAGAGAAGTCTGTACATATGTCAACACTGAGGGTAGAGCTCAGCAGACTAAGGTAGCAGTGACA
CCTCCTGGCTTGCAAGAGAAGACTGGAAAATTATAAGAGCACTCTCTGAGATTGCTGGAATGACTCTT
CCATATGATACTTGATCAAGTAAGGAACAGATTGGAAGAAGTCTCTCCTAATCTTGTTCGATATGAT
GATATTGAAGGGGCTAATTACTTCCAGCAAGCAAATGAGCTCTCAAAGCTAGTGAACCAGCAGCTTCTT
GCTGACCACTTGTTCACCTCAGCTAACTATAAAAGACTTCTACATGACAGATTCAATTAGCAGAGCC
TCACAGACAATGGCCAAATGTGTCAAAGCTGTACAGAGGGTGCCAGGCAGTAGAGGAACCATCCATA
TGCTGA
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_001199981

Insert Size: 2076 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001199981.1](#)

RefSeq Size: 3365 bp

RefSeq ORF: 2076 bp

Locus ID: 4719

UniProt ID: [P28331](#)

Cytogenetics: 2q33.3

Protein Pathways: Alzheimer's disease, Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease

MW: 75.4 kDa

Gene Summary: The protein encoded by this gene belongs to the complex I 75 kDa subunit family. Mammalian complex I is composed of 45 different subunits. It locates at the mitochondrial inner membrane. This protein has NADH dehydrogenase activity and oxidoreductase activity. It transfers electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone. This protein is the largest subunit of complex I and it is a component of the iron-sulfur (IP) fragment of the enzyme. It may form part of the active site crevice where NADH is oxidized. Mutations in this gene are associated with complex I deficiency. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2011]

Transcript Variant: This variant (2) lacks an alternate in-frame exon compared to variant 1. The resulting isoform (2) has the same N- and C-termini but is shorter compared to isoform 1.