

## Product datasheet for **SC329675**

### **CUX1 (NM\_001202546) Human Untagged Clone**

#### **Product data:**

<b>Product Type:</b>	Expression Plasmids
<b>Product Name:</b>	CUX1 (NM_001202546) Human Untagged Clone
<b>Tag:</b>	Tag Free
<b>Symbol:</b>	CUX1
<b>Synonyms:</b>	CASP; CDP; CDP/Cut; CDP1; Clox; COY1; CUTL1; CUX; Cux/CDP; GDDI; GOLIM6; Nbla10317; p75; p100; p110; p200
<b>Vector:</b>	pCMV6-Entry (PS100001)



[View online »](#)

**Fully Sequenced ORF:** >SC329675 representing NM\_001202546.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

ATGGCGGCCAATGTGGGATCGATGTTTCAATATTGGAAGCGCTTTGATTTACAGCAGCTGCAGGATTTG
CGAAGCAGGTAGCGCCGCTGCTGAAGAGTTTCCAAGGAGAGATTGATGCACTGAGTAAAAGAAGCAAG
GAAGCTGAAGCAGCTTTCTTGAATGTCTACAAAAGATTGATTGACGTCCCAGATCCCCTACCAGCTTTG
GATCTCGGACAGCAACTCCAGCTCAAAGTGCAGCGCTGCACGATATTGAAACAGAGAACCAGAACTT
AGGAAAACCTCTGGAAGAATAACAACAAGGAATTTGCTGAAGTGAAAAATCAAGAGGTTACGATAAAAGCA
CTTAAAGAGAAAAATCCGAGAATATGAACAGACACTGAAGAACCAAGCCGAAACCATAGCTCTTGAGAAG
GAACAGAAGTTACAGAATGACTTTGCAGAAAAGGAGAGAAAGCTGCAGGAGACACAGATGTCCACCACC
TCAAAGCTGGAGGAAGCTGAGCATAAGGTTTCAGAGCCTACAAACAGCCCTGGAAAAAACTCGAACAGAA
TTATTTGACCTGAAAACCAAATACGATGAAGAACTACTGCAAAGCCGACGAGATTGAAATGATCATG
ACGGACCTTGAAGGGCAAACCAGAGGGCAGAGTGGCTCAGAGAGAGGGGAGACCTTAAGGGAACAG
CTCTCATCGGCCAATCACTCCCTCCAGCTGGCCTCACAGATCCAGAAGGCCAGACGTGGCCATAGAG
GTGCTGACCCGCTCCAGCCTAGAAGTTGAGTTGGCCGCAAGGAGCGGGAGATCGCACAGCTGGTGGAG
GACGTGCAGAGACTCCAGGCCAGCCTACCAAGCTGCGGGAGAATTTCGGCCAGCCAGATCTCACAGCTT
GAGCAGCAGCTGAGCGCCAAAAACAGCACACTCAAACAACCTGGAAGAAAACTCAAAGGCCAGGCTGAC
TATGAAGAGGTGAAGAAAGAGCTGAACATTCTGAAGTCCATGGAGTTTGCACCGTCCGAGGGCGCTGGG
ACACAGGATGCGGCCAAGCCCCTGGAGGTGCTGTTGCTGGAGAAGAACCCTCGTGCAGTCCGAGAAC
GCCGCGCTGCGCATCTCCAACAGCGACCTGAGCGGACGCTGTGCAGAGCTGCAAGTCCGTATCACTGAG
GCTGTGGCCACAGCCACTGAGCAGAGAGAGCTGATCGCCCGCTGGAGCAGGACCTGAGCATCATTAG
TCCATCCAGCGGCCGATGCCGAGGGTCCGCTGAGCACCCTGGAGAAGATCCAGAGCCCATCAA
GAGGCCACTGCCCTATTCTACGACCTGCAGCACCAGCCAGCGGTGCCCTCCAGAGGGCCAGTGGAT
TCACTGCTTTCCATCATCTCCAGCCAGAGGGAGCGCTTCCGTGCCGGAACCAGGAGCTTGAGGCCGAG
AACCGCCTGGCCAGCACACCCTCCAGGCCCTGCAGAGTGAAGTGGACAGCCTGCGCGCCGACAACATC
AAGCTCTTTGAGAAGATCAAGTTCCTGCAGAGCTACCCTGGCCGGGGCAGCGGCAGTGATGACACGGAG
CTGCGGTACTCGTCCCAGTACGAGGAGCGCTGGACCCCTTCTCCTCCTCAGCAAGCGGGAGCGGCAG
AGGAAGTACCTGAGCTTGAGTCCCTGGGACAAGGCCACCCTCAGCATGGGGCGTCTGTTCTCTCCAAC
AAGATGGCGCGCACCATCGGCTTCTTACACACTGTTCTGCACTGCCTGGTCTTCTGGTGTCTAC
AAGCTGGCATGGAGCGAGAGCATGGAGAGGGACTGTGCCACCTTCTGCGCCAAGAAGTTCGCTGACCAC
CTGCACAAGTTCACGAGAATGACAACGGGGCTGCGGCTGGTACTTGTGGCAGTGA
  
```

**Restriction Sites:**

Sgfl-Mlul

**Plasmid Map:**



<b>ACCN:</b>	NM_001202546
<b>Insert Size:</b>	1920 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001202546.1</a>
<b>RefSeq Size:</b>	2818 bp
<b>RefSeq ORF:</b>	1920 bp
<b>Locus ID:</b>	1523
<b>UniProt ID:</b>	<a href="#">Q13948</a>
<b>Cytogenetics:</b>	7q22.1
<b>Protein Families:</b>	Transcription Factors, Transmembrane
<b>MW:</b>	72.8 kDa
<b>Gene Summary:</b>	<p>The protein encoded by this gene is a member of the homeodomain family of DNA binding proteins. It may regulate gene expression, morphogenesis, and differentiation and it may also play a role in the cell cycle progression. Several alternatively spliced transcript variants encoding different isoforms have been identified.[provided by RefSeq, Feb 2011]</p> <p>Transcript Variant: This variant (7) lacks an in-frame exon in the 5' region and has an alternate 3' sequence including the coding region, as compared to variant 1. The resulting isoform (g) lacks an internal segment and has a shorter and distinct C-terminus, as compared to isoform d.</p>