

Product datasheet for **SC329681**

Ribosomal protein S10 (RPS10) (NM_001203245) Human Untagged Clone


Product data:

Product Type:	Expression Plasmids
Product Name:	Ribosomal protein S10 (RPS10) (NM_001203245) Human Untagged Clone
Tag:	Tag Free
Symbol:	RPS10
Synonyms:	DBA9; S10
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC329681 representing NM_001203245. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGTTGATGCCTAAGAAGAACCGGATTGCCATTTATGAACTCCTTTTTAAGGAGGGAGTCATGGTGGCC
AAGAAGGATGTCCACATGCCTAAGCACCCGGAGCTGGCAGACAAGAATGTGCCAACCTTCATGTCATG
AAGGCCATGCAGTCTCTCAAGTCCCAGGCTACGTGAAGGAACAGTTTGCCTGGAGACATTTCTACTGG
TACCTTACCAATGAGGGTATCCAGTATCTCCGTGATTACCTTCATCTGCCCCGGAGATTGTGCCTGCC
ACCCTACGCCGTAGCCGTCCAGAGACTGGCAGGCCTCGGCCTAAAGGTCTGGAGGGTGAGCGACCTGCC
AGACTCACAAGAGGGGAAGCTGACAGAGATACCTACAGACGGAGTGCTGTGCCACCTGGTGCCGACAAG
AAAGCCGAGGCTGGGGCTGGGTGAGCAACCGAATTCCAGTTTAGAGGCGGATTTGGTCGTGGACGTGGT
CAGCCACCTCAGTAA
```

Restriction Sites: SgfI-MluI

Plasmid Map:



[View online »](#)

ACCN:	NM_001203245
Insert Size:	498 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001203245.2
RefSeq Size:	813 bp
RefSeq ORF:	498 bp
Locus ID:	6204
UniProt ID:	P46783
Cytogenetics:	6p21.31
Protein Pathways:	Ribosome
MW:	18.9 kDa

Gene Summary:

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 40S subunit. The protein belongs to the S10E family of ribosomal proteins. It is located in the cytoplasm. Variable expression of this gene in colorectal cancers compared to adjacent normal tissues has been observed, although no correlation between the level of expression and the severity of the disease has been found. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. Alternate splicing results in multiple transcript variants that encode the same protein. Naturally occurring read-through transcription occurs between this locus and the neighboring locus NUDT3 (nudix (nucleoside diphosphate linked moiety X)-type motif 3).[provided by RefSeq, Feb 2011]

Transcript Variant: This variant (1) represents the longer transcript. Variants 1, 2 and 3 encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.