

Product datasheet for **SC329701**

ANKS1B (NM_001204081) Human Untagged Clone

Product data:

| | |
|----------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | ANKS1B (NM_001204081) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | ANKS1B |
| Synonyms: | AIDA; AIDA-1; ANKS2; cajalin-2; EB-1; EB1 |
| Vector: | pCMV6-Entry (PS100001) |
| Fully Sequenced ORF: | >SC329701 representing NM_001204081. Blue=Insert sequence Red=Cloning site Green=Tag(s) |

```
ATGCAGGGCGATGCAAGGAGGAGAAGAAATGAAACTACTTTGATGATATCCCCGATCAAACTGGAG
AGGCAGATGGCTCAGACAGGAGACTGGGGAGAACCTTCCATTACCTTGCGACCTCCGAATGAAGCCACA
GCCTCTACCCCGGTACAGTACTGGCAGCATCACCCAGAAAAGCTTATCTTCCAGTCGTGTGATTACAAA
GCTTTTTATTTAGGTTCTATGCTGATAAAAGAGCTTAGGGGGACAGAATCAACCCAAGATGCTTGTGCA
AAAATGCGGGCTAACTGTCAGAAGTCTACAGAGCAAATGAAGAAGTCCCTACTATTATTCTTTCTGTC
TCATATAAAGGAGTCAAATTTATTGATGCAACAAATAAGAACATAAATTGCTGAGCATGAAATTCGTAAT
ATCTCCTGTGCTGCCAGGACCCAGAAGACCTCTCAACATTTGCCTATATCACAAAAGATTTGAAGTCT
AATCACCCTACTGTCATGTGTTTACTGCCTTTGATGTGAATTTAGCCTATGAAATCATCCTAACCCCTG
GGACAGGCATTCGAAGTCGTTACCAGCTAGCACTACAAGCAAGAAAAGGGGGACTCCTCCACTT
CCAGAAAGCTTTGAAAACAAACCCTCAAACCCATCCCCAAGCCCCGCGTTAGCATTGCAAGTCCGTG
GATCTTCTTCATGCCAGCCATACCGGGCAGGAACCATCAGAGAGACACTGAGGAGGCACTCAGAAAA
TTTTGA
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:


ACCN: NM_001204081

Insert Size: 765 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001204081.1](#)

RefSeq Size: 2043 bp

RefSeq ORF: 765 bp

Locus ID: 56899

UniProt ID: [Q7Z6G8](#)

Cytogenetics: 12q23.1

MW: 28.9 kDa

Gene Summary: This gene encodes a multi-domain protein that is predominantly expressed in brain and testis. This protein interacts with amyloid beta protein precursor (AbetaPP) and may have a role in normal brain development, and in the pathogenesis of Alzheimer's disease. Expression of this gene has been shown to be elevated in patients with pre-B cell acute lymphocytic leukemia associated with t(1;19) translocation. Alternatively spliced transcript variants encoding different isoforms (some with different subcellular localization, PMID:15004329) have been described for this gene. [provided by RefSeq, Aug 2011]
Transcript Variant: This variant (12) differs in the 5' UTR and coding region compared to variant 1. The resulting isoform (l) has a shorter and distinct N-terminus compared to isoform a.