

Product datasheet for **SC329717**

SGK3 (C8orf44-SGK3) (NM_001204173) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: SGK3 (C8orf44-SGK3) (NM_001204173) Human Untagged Clone
Tag: Tag Free
Symbol: C8orf44-SGK3
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC329717 representing NM_001204173.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

ATGCAAAGAGATCACACCATGGACTACAAGGAAAGCTGCCAAGTGAAGCATTCCCAGCTCCGATGAA
CACAGAGAGAAAAAGAAGAGGTTTACTGTTTATAAAGTTCTGGTTTCAGTGGGAAGAAGTGAATGGTTT
GTCTTCAGGAGATATGCAGAGTTTGATAAACTTTATAACACTTTAAAAAACAGTTTCTGCTATGGCC
CTGAAGATTCTGCCAAGAGAATATTTGGTGATAATTTTGATCCAGATTTTATTAACAAGACGAGCA
GGACTAAACGAATTCATTCAGAACCTAGTTAGGTATCCAGAACTTTATAACCATCCAGATGTCAGAGCA
TTCCCTTCAAATGGACAGTCCAAAACACCAGTCAGATCCATCTGAAGATGAGGATGAAAGAAGTTCTCAG
AAGCTACACTCTACCTCACAGAACATCAACCTGGGACCGTCTGGAAATCCTCATGCCAAACCAACTGAC
TTTGATTTCTTAAAAGTTATTGGAAAAGCAGCTTTGGCAAGGTTCTTCTTGCAAAACGGAACTGGAT
GGAAAATTTTATGCTGTCAAAGTGTACAGAAAAAATAGTTCTCAACAGAAAAGAGCAAAAACATATT
ATGGCTGAACGTAATGTGCTCTTGAAAAATGTGAAACATCCGTTTTTGGTTGGATTGCATTATTCCTTC
CAAACAACTGAAAAGCTTTATTTTGTCTGGATTTTGTTAATGGAGGGGAGCTTTTTTCCACTTACAA
AGAGAACGGTCTTTCTGAGCACAGAGCTAGGTTTTACGCTGCTGAAATGCTAGTGCATTGGGTAC
TTACATTCATCAAAATAGTATACAGAGACTTGAAACCAGAAAATATTCTTTTGGATTGAGTAGGACAT
GTTGTCTTAACAGATTTTGGGCTTTGTAAGAAGGAATTGCTATTTCTGACACCACTACCACATTTGT
GGGACACCAGAGTATCTTGACCTGAAGTAATTAGAAAACAGCCCTATGACAATACTGTAGATTGGTGG
TGCTTGGGGCTGTTCTGTATGAAATGCTGTATGGATTGCCTCCTTTTTATTGCCGAGATGTTGCTGAA
ATGTATGACAATATCCTTCACAAACCCCTAAGTTTGAGGCCAGGAGTGAGTCTTACAGCCTGGTCCATT
CTGGAAGAAGTCTAGAAAAAGACAGGCAAAATCGACTTGGTGCCAAGGAAGACTTTCTTGAAATTCAG
AATCATCCTTTTTTTGAATCACTCAGCTGGGCTGACCTTGTACAAAAGAAGATCCACCACCATTTAAT
CCTAATGTGGCTGGACCAGATGATATCAGAACTTTGACACAGCATTTACAGAAGAAACAGTTCCATAT
TCTGTGTGTGATCTTCTGACTATTCTATAGTGAATGCCAGTGTATTGGAGGCAGATGATGCATTCGTT
GGTTTCTTATGCACCTCCTTCAGAAGACTTATTTTTGTGA
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001204173

Insert Size: 1491 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001204173.1](#)

RefSeq Size: 4563 bp

RefSeq ORF: 1491 bp

Locus ID: 100533105

UniProt ID: [Q96BR1](#)

Cytogenetics: 8q13.1

MW: 57.1 kDa

Gene Summary: This locus represents naturally occurring read-through transcription between the neighboring putative uncharacterized protein C8orf44 (GeneID 56260) and serine/threonine-protein kinase Sgk3 (GeneID 23678) genes on chromosome 8. The read-through transcript produces a protein that shares sequence identity with the downstream gene product. [provided by RefSeq, Feb 2011]