

## Product datasheet for **SC329744**

### **SAP97 (DLG1) (NM\_001204386) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	SAP97 (DLG1) (NM_001204386) Human Untagged Clone
Tag:	Tag Free
Symbol:	DLG1
Synonyms:	DLGH1; hdlg; SAP-97; SAP97
Vector:	pCMV6-Entry (PS100001)




[View online »](#)

Fully Sequenced ORF: >SC329744 representing NM\_001204386.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGCCGGTCCGGAAGCAAGATACCCAGAGAGCATTGCACCTTTTGGAGGAATATCGTTCAAACCTAAGC
CAAACCTGAAGACAGACAGCTCAGAAGTTCATAGAACGGGTTATTAACATATTTCCAGAGCAACCTCTTT
CAGGCTTTAATAGATATCAAGAATTTTATGAAGTGACCTTACTGGATAATCCAAATGTATAGATCGT
TCAAAGCCGTCTGAACCAATCAACCTGTGAATACTTGGGAGATTTCCAGCCTTCCAAGCTACTACTGTG
ACTTCAGAGACACTGCCAAGCAGCCTTAGCCCTAGTGTAGAGAAATACAGGTATCAGGATGAAGATACA
CCTCCTCAAGAGCATATTTCCCCACAAATCACAAATGAAGTGATAGGTCCAGAATTGGTTCATGTCTCA
GAGAAGAACTTATCAGAGATTGAGAATGTCCATGGATTTGTTTCTCATTCTCATATTTCCACCAATAAAG
GCAAACTCTCCCCAGTACTGGTCAACACAGATAGCTTGGAAACACCAACTTACGTTAATGGCACAGAT
GCAGATTATGAATATGAAGAAATCACACTTGAAAGGGGAAATTCAGGGCTTGGTTTCAGCATTGCAGGA
GGTACGGACAACCCACACATTGGAGATGACTCAAGTATTTTCATTACCAAAATTATCACAGGGGGAGCA
GCCGCCAAGATGGAAGATTGCGGGTCAATGACTGTATATTACGAGTAAATGAAGTAGATGTTTCGTGAT
GTAACACATAGCAAAGCAGTTGAAGCGTTGAAAGAAGCAGGGTCTATTGTACGCTTGTATGTA AAAAGA
AGGAAACCAAGTGTGCAAAAAAATATGGAATAAAGCTCATTAAAGGTCCATAAAGTCTTGGGTTTAGC
ATTGCTGGAGGTGTTGGAATCAGCATATTCCTGGGGATAATAGCATCTATGTAACCAAAAATAATTGAA
GGAGGTGCAGCACATAAGGATGGCAAACTTCAGATTGGAGATAAACTTTTAGCAGTGAATAACGTATGT
TTAGAAGAAGTACTCATGAAGAAGCAGTAAGTGCCTTAAAGAACACATCTGATTTTGTATTTGAAA
GTGGCAAACCCACAAGTATGTATATGAATGATGGCTATGCACCACCTGATATCACCAACTCTTCTTCT
CAGCCTGTTGATAACCATGTTAGCCCATCTTCTTCTTGGGCCAGACACCAGCATCTCCAGCCAGATAC
TCCCAGTTTCTAAAGCAGTACTTGGAGATGATGAAATTAACAAGGAACTAGAAAAGTGTGTTCTTCACT
CGTGGCTCAACGGGCCTTGGTTTCAACATTGTAGGAGGAGAAGATGGAGAAGGAATATTTATTTCTTTT
ATCTTAGCCGGAGGACCTGCTGATCTAAGTGGAGAGCTCAGAAAAGGAGATCGTATTATATCGGTA AAC
AGTGTTGACCTCAGAGCTGCTAGTCATGAGCAGGAGCAGCTGCATTGAAAAATGCTGGCCAGGCTGTC
ACAATTGTTGCACAATATCGACCTGAAGAATACAGTCGTTTTGAAAGCTAAAATACATGATTTACGGGAG
CAGATGATGAATAGTAGTATTAGTTCAGGGTCAGGTTCTCTTTCGAACTAGCCAGAAGCGATCCCTCTAT
GTCAGAGCCCTTTTTGATTATGACAAGACTAAAGACAGTGGGCTTCCCAGTCAGGGACTGAACTTCAA
TTTGGAGATATCCTCCATGTTATTAATGCTTCTGATGATGAATGGTGGCAAGCCAGGCAGGTTACACCA
GATGGTGAGAGCGATGAGGTCGGAGTATTCCCAGTAAACGCAGAGTTGAGAAGAAAGAACGAGCCCGA
TTAAAAACAGTGAATTC AATTCTAAAACGAGAGATAAAGGGTCATTCAATGACAAGCGTAAAAAGAAC
CTTTTTCCCGAAAATCCCCTTCTACAAGAACAAGGACCAGAGTGAGCAGGAAACAAGTATGCTGAC
CAGCATGTAACCTTAATGCCAGCGATAGTGAAAGTAGTTACCGTGGTCAAGAAGATACGCTTATCT
TATGAACCAAGTGAATCAACAAGAAGTAAATTATACTCGACCAGTGATCATATTGGGACCTATGAAAGAC
AGGATAAATGATGACTTGATCTCAGAAATTTCTGACAAATTTGGATCCTGTGTTCCCTCATACAAC TAGA
CCAAAACGAGATTATGAGGTAGATGGAAGAGATTATCATTTTGTGACTTCAAGAGAGCAGATGGAAAA
GATATCCAGGAACATAAATTCATTGAAGCTGGCCAGTATAACAATCATCTATATGGAACAAGTGTTCAG
TCTGTACGAGAAGTAGCAGAAAAGGGCAACACTGTATCCTTGATGTGTCTGGAATGCCATAAAGAGA
TTACAGATTGCACAGCTTACCCTATCTCCATTTTATTAACCCAAATCCATGGAAAATATCATGGAA
ATGAATAAGCGTCTAACAGAAGAACAAGCCAGAAAAACATTTGAGAGAGCCATGAAACTGGAACAGGAG
TTTACTGAACATTTACAGCTATTGTACAGGGGATACGCTGGAAGACATTTACAACCAAGTGAACAG
ATCATAGAAGAACAATCTGGTCTTACATCTGGGTTCCGGCAAAGAAAAGCTATGA
```

Restriction Sites: Sgfl-Mlul

**Plasmid Map:**


**ACCN:** NM\_001204386

**Insert Size:** 2679 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001204386.1](#)

**RefSeq Size:** 4941 bp

**RefSeq ORF:** 2679 bp

**Locus ID:** 1739

**UniProt ID:** [Q12959](#)

**Cytogenetics:** 3q29

**Protein Families:** Druggable Genome

**Protein Pathways:** T cell receptor signaling pathway

**MW:** 99.8 kDa

**Gene Summary:** This gene encodes a multi-domain scaffolding protein that is required for normal development. This protein may have a role in septate junction formation, signal transduction, cell proliferation, synaptogenesis and lymphocyte activation. A multitude of transcript variants deriving from alternative splicing and the use of multiple alternate promoter have been observed, including some splice variants that may be specific to brain and other tissues. An upstream uORF may regulate translation at some splice variants of this gene. [provided by RefSeq, Sep 2018]

Transcript Variant: This variant (3) is missing an in-frame coding exon, and contains an alternate in-frame coding exon compared to variant 1. This results in a shorter isoform (3) compared to isoform 1.