

Product datasheet for **SC329746**

SAP97 (DLG1) (NM_001204388) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SAP97 (DLG1) (NM_001204388) Human Untagged Clone
Tag:	Tag Free
Symbol:	DLG1
Synonyms:	DLGH1; hdlg; SAP-97; SAP97
Vector:	pCMV6-Entry (PS100001)




[View online »](#)

Fully Sequenced ORF: >SC329746 representing NM_001204388.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGAATTACATTTTCGGAAACAACACACTTCTATACTCTCGCGGCAGTCGAGGAGGCAATACTAGCTCT
AGCCATGGCTCAGCAGGCCAAAGCAGAAACACTGGGCAAAAAGGGCTCGTCAGATGAACTGCAAGCT
GAGCCAGAACCTTCACGCTGGCAGCAGATAGTTGCATTTTTCACTCGAAGACACAGCTTTATTGACTGC
ATCTCGGTAGCCACCAGCTCCACCAGGCAAATCCTCCCCAGTACTGGTCAACACAGATAGCTTGGAA
ACACCAACTTACGTTAATGGCACAGATGCAGATTATGAATATGAAGAAATCACACTTGAAAGGGGAAAT
TCAGGGCTTGGTTTCAGCATTGCAGGAGGTACGGACAACCCACACATTGGAGATGACTCAAGTATTTTC
ATTACAAAATTATCACAGGGGGAGCAGCCGCCAAGTGAAGATTGCGGGTCAATGACTGTATATTA
CGAGTAAATGAAGTAGATGTTCTGTGATGAACACATAGCAAAGCAGTTGAAGCGTTGAAAGAAGCAGGG
TCTATTGTACGCTTGTATGTAAGAAAGAAAGGAAACCAGTGTGAGAAAAAATATGGAATAAAGCTCATT
AAAGGTCCTAAAGGCTTGGGTTTAGCATTGCTGGAGGTGTTGGAATCAGCATATTCCTGGGGATAAT
AGCATCTATGTAACCAAAATAATTGAAGGAGGTGCAGCACATAAGGATGGCAAACCTCAGATTGGAGAT
AACTTTTAGCAGTGAATAACGTATGTTAGAAGAAGTACTCATGAAGAAGCAGTAACTGCCTTAAG
AACACATCTGATTTTGTATTTGAAAGTGGCAAAACCCACAAGTATGTATATGAATGATGGCTATGCA
CCACCTGATATCACCAACTCTTCTCTCAGCCTGTTGATAACCATGTTAGCCCATCTTCTTCTTGGGC
CAGACACCAGCATCTCCAGCCAGATACTCCCCAGTTTCTAAAGCAGTACTTGGAGATGATGAAATTACA
AGGGAACCTAGAAAAGTTGTTCTTCATCGTGGCTCAACGGGCCTTGGTTTCAACATTGTAGGAGGAGAA
GATGGAGAAGGAATATTTATTTCTTTATCTTAGCCGGAGGACCTGCTGATCTAAGTGGAGAGCTCAGA
AAAGGAGATCGTATTATATCGGTAACAGTGTGACCTCAGAGCTGCTAGTCATGAGCAGGCAGCAGCT
GCATTGAAAAATGCTGGCCAGGCTGTACAATTGTTGCACAATATCGACCTGAAGAATACAGTCGTTTT
GAAGCTAAAATACATGATTTACGGGAGCAGATGATGAATAGTAGTATTAGTTCAAGGTCAGGTTCTCTT
CGAACTAGCCAGAAGCGATCCCTCTATGTGAGAGCCCTTTTTGATTATGACAAGACTAAAGACAGTGGG
CTTCCCAGTCAGGGACTGAACTTCAAATTTGGAGATATCCTCCATGTTATTAATGCTTCTGATGATGAA
TGGTGGCAAGCCAGGCAGGTTACACCAGATGGTGAGAGCGATGAGGTGGGAGTATTCCCAGTAAACGC
AGAGTTGAGAAGAAAGAACGAGCCCGATTA AAAACAGTGAATTC AATTCTAAAACGAGAGATAAAGGG
GAGATCCCTGACGACATGGGATCAAAGGCCTGAAGCATGTAACCTCTAATGCCAGCGATAGTAAAAGT
AGTTACCGTGGTCAAGAAGAATACGTCTTATCTTATGAACCAGTGAATCAACAAGAAGTTAATTACT
CGACCAGTGATCATATTGGGACCTATGAAAGACAGGATAAATGATGACTTGATCTCAGAATTTCTGAC
AAATTTGGATCCTGTGTTCTCATACAACTAGACCAAAACGAGATTATGAGGTAGATGGAAGAGATTAT
CATTTTGTGACTTCAAGAGAGCAGATGGAAAAAGATATCCAGGAACATAAATTCATTGAAGCTGGCCAG
TATAACAATCATCTATATGGAACAAGTGTTCAGTCTGTACGAGAAGTAGCAGAAAAGGGCAAACACTGT
ATCCTTGATGTGCTGGAAATGCCATAAAGAGATTACAGATTGCACAGCTTTACCCTATCTCCATTTTT
ATTAACCCAAATCCATGGAAAATATCATGGAAATGAATAAGCGTCTAACAGAAGAACAAGCCAGAAAA
ACATTTGAGAGAGCCATGAAACTGGAACAGGAGTTTACTGAACATTTACAGCTATTGTACAGGGGGAT
ACGCTGGAAGACATTTACAACCAAGTGAACAGATCATAGAAGAACAATCTGGTTCTTACATCTGGGT
CCGGCAAAAGAAAAGCTATGA
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001204388

Insert Size: 2367 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001204388.1](#)

RefSeq Size: 4762 bp

RefSeq ORF:	2367 bp
Locus ID:	1739
UniProt ID:	Q12959
Cytogenetics:	3q29
Protein Families:	Druggable Genome
Protein Pathways:	T cell receptor signaling pathway
MW:	87.2 kDa
Gene Summary:	<p>This gene encodes a multi-domain scaffolding protein that is required for normal development. This protein may have a role in septate junction formation, signal transduction, cell proliferation, synaptogenesis and lymphocyte activation. A multitude of transcript variants deriving from alternative splicing and the use of multiple alternate promoter have been observed, including some splice variants that may be specific to brain and other tissues. An upstream uORF may regulate translation at some splice variants of this gene. [provided by RefSeq, Sep 2018]</p> <p>Transcript Variant: This variant (5) contains an alternate 5' terminal exon and is missing an in-frame coding exon compared to variant 1. This results in a shorter isoform (5) with a unique N-terminus compared to isoform 1.</p>