

## Product datasheet for **SC329750**

### **XIAP (NM\_001204401) Human Untagged Clone**

#### Product data:

Product Type:	Expression Plasmids
Product Name:	XIAP (NM_001204401) Human Untagged Clone
Tag:	Tag Free
Symbol:	XIAP
Synonyms:	API3; BIRC4; hIAP-3; hIAP3; IAP-3; ILP1; MIHA; XLP2
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC329750 representing NM_001204401. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```


ATGACTTTTAACAGTTTTGAAGGATCTAAAACCTGTGTACCTGCAGACATCAATAAGGAAGAAGAATTT
GTAGAAGAGTTTAAATAGATTA AAAA ACTTTTGCTAATTTTCCAAGTGGTAGTCCTGTTTCAGCATCAACA
CTGGCAGCAGCAGGGTTTCTTTACTGGTGAAGGAGATACCGTGC GG TGCTTTAGTTGTCATGCAGCT
GTAGATAGATGGCAATATGGAGACTCAGCAGTTGGAAGACACAGGAAAGTATCCCCAAATGCAGATTT
ATCAACGGCTTTTATCTTGAAAATAGTGCCACGCAGTCTACAAAATCTGGTATCCAGAATGGTCAGTAC
AAAGTTGAAAACCTCTGGGAAGCAGAGATCATTTTGCCTTAGACAGGCCATCTGAGACACATGCAGAC
TATCTTTTGAGAACTGGGCAGGTTGTAGATATATCAGACACCATATACCCGAGGAACCCCTGCCATGTAT
AGTGAAGAAGCTAGATTAAGATCCTTTCAGAACTGGCCAGACTATGCTCACCTAACCCCAAGAGAGTTA
GCAAGTGCTGGACTCTACTACACAGGTATTGGTGACCAAGTGCAGTGCTTTTGTGTGGTGGAAAACCTG
AAAAATGGGAACCTTGTGATCGTGCCTGGTCAGAACACAGGCGACACTTTCCTAATTGCTTCTTTGTT
TTGGGCCGGAATCTTAATATTGGAAGTGAATCTGATGCTGTGAGTTCTGATAGGAATTTCCCAAATTC
ACAAATCTTCCAAGAAATCCATCCATGGCAGATTATGAAGCACGGATCTTTACTTTTGGGACATGGATA
TACTCAGTTAACAAGGAGCAGCTTCAAGAGCTGGATTTTATGCTTTAGGTGAAGGTGATAAAGTAAAG
TGCTTTCACTGTGGAGGAGGGCTAACTGATTGGAAGCCAGTGAAGACCCTTGGGAACAACATGCTAAA
TGGTATCCAGGGTGCAAATATCTGTTAGAACAAGAGGACAAGAATATATAAACAATATTCAATTTAACT
CATTCACTTGAGGAGTGTCTGGTAAGAATACTGAGAAAACCATCACTAACTAGAAGAATTGATGAT
ACCATCTTCCAAAATCCTATGGTACAAGAAGCTATACGAATGGGGTTCAGTTTCAAGGACATTAAGAAA
ATAATGGAGGAAAAAATTCAGATATCTGGGAGCAACTATAAATCACTTGAGGTTCTGTTGCAGATCTA
GTGAATGCTCAGAAAGACAGTATGCAAGATGAGTCAAGTCAGACTTCATTACAGAAAGAGATTAGTACT
GAAGAGCAGCTAAGGCGCCTGCAAGAGGAGAAGCTTTGCAAAAATCTGTATGGATAGAAAATTTGCTATC
GTTTTTGTTCCTTGTGGACATCTAGTCACTTGTAAACAATGTGCTGAAGCAGTTGACAAGTGTCCCATG
TGCTACACAGTCATTACTTTCAAGCAAAAAATTTTATGTCTAA

```

Restriction Sites: SgfI-MluI



[View online »](#)

**Plasmid Map:**


**ACCN:** NM\_001204401

**Insert Size:** 1494 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001204401.1](#)

**RefSeq Size:** 8427 bp

**RefSeq ORF:** 1494 bp

**Locus ID:** 331

**UniProt ID:** [P98170](#)

<b>Cytogenetics:</b>	Xq25
<b>Protein Families:</b>	Druggable Genome
<b>Protein Pathways:</b>	Apoptosis, Focal adhesion, NOD-like receptor signaling pathway, Pathways in cancer, Small cell lung cancer, Ubiquitin mediated proteolysis
<b>MW:</b>	56.7 kDa
<b>Gene Summary:</b>	<p>This gene encodes a protein that belongs to a family of apoptotic suppressor proteins. Members of this family share a conserved motif termed, baculovirus IAP repeat, which is necessary for their anti-apoptotic function. This protein functions through binding to tumor necrosis factor receptor-associated factors TRAF1 and TRAF2 and inhibits apoptosis induced by menadione, a potent inducer of free radicals, and interleukin 1-beta converting enzyme. This protein also inhibits at least two members of the caspase family of cell-death proteases, caspase-3 and caspase-7. Mutations in this gene are the cause of X-linked lymphoproliferative syndrome. Alternate splicing results in multiple transcript variants. Pseudogenes of this gene are found on chromosomes 2 and 11.[provided by RefSeq, Feb 2011]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Both variants 1 and 2 encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>