

Product datasheet for **SC329790**

CRISP1 (NM_001205220) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: CRISP1 (NM_001205220) Human Untagged Clone
Tag: Tag Free
Symbol: CRISP1
Synonyms: AEGL1; ARP; CRISP-1; HEL-S-57; HSCRISP1D; HSCRISP1G; HUMARP
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC329790 representing NM_001205220.
Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
ATGGAATTAACACCTCTGTTTTGGTTGCTGCTGCTTACTGCCTATGTTGCCATGAAAAAG
AAATCAGCTAGAGACCAATTTAATAAGCTCGTCACCGACTTGCCAAATGTACAAGAAGAGATCGTTAAT
ATACACAACGCCCTCAGGAGAAGAGTAGTTCCACCAGCCAGCAACATGCTGAAGATGAGTTGGAGTGAA
GAGGCTGCACAAAATGCCAGAATTTTTCAAAGTATTGTGATATGACAGAGAGCAACCCCTTGAGAGG
AGACTTCAAATACCTTTTGGGAGAAAATATGCATATGACATCTTATCCTGTATCATGGTCAAGTGTA
ATTGGAGTCTGGTACAGTGAGTCTACAAGTTTCAAACATGGAGAATGGACAACAACGGATGATGACATA
ACTACTGACCACTACACTCAGATTGTTTGGGCCACATCTTACCTGATTGGCTGTGCCATTGCATCTTGC
CGCCAACAAGGATCACCTCGATATCTCTACGTTTGTCACTATTGTGATGAGGAAATGATCCTGAAACA
AAGAATGAACCTTATAAGACAGGCGTCCCATGTGAAGCCTGCCCAAGTAAGTGTGAAGACAACTTTGC
ACTAACCCCTGCATCTACTATGATGAATACTTCGACTGTGACATACAAGTCCATTATCTGGGATGCAAC
CACTCAACAACATCCTATTCTGTAAAGCCACTTGCTGTGTGACTGAGATAAAATAG
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:



ACCN: NM_001205220

Insert Size: 750 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001205220.1](#)

RefSeq Size: 1903 bp

RefSeq ORF: 750 bp

Locus ID: 167

UniProt ID: [P54107](#)

Cytogenetics: 6p12.3

MW: 28.5 kDa

Gene Summary: Fertilization consists of a sequence of specific cell-cell interactions culminating in the fusion of the sperm and egg plasma membranes. Recognition, binding, and fusion occur through the interaction of complementary molecules that are localized to specific domains of the sperm and egg plasma membranes. In the sperm, the postacrosomal region or equatorial segment is involved in sperm-egg plasma membrane fusion. The protein encoded by this gene is a member of the cysteine-rich secretory protein (CRISP) family. It is expressed in the epididymis, is secreted into the epididymal lumen, and binds to the postacrosomal region of the sperm head, where it plays a role in sperm-egg fusion. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2011]

Transcript Variant: This variant (3) contains an alternate 5' terminal non-coding exon, thus differs in the 5' UTR compared to variant 1. Variants 1 and 3 encode the same isoform (1).