

Product datasheet for **SC329799**

THEG5 (NM_001205273) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: THEG5 (NM_001205273) Human Untagged Clone
Tag: Tag Free
Symbol: THEG5
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC329799 representing NM_001205273.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

ATGGCCATGCATTTATCTGGAGTCAATGCTGATTCTAACTCCATGGCCACCTTTACCCTCAGCTTCATT
 GTAGCTGTGAATGCTCCGTTCTCTCTTTTACCCTGAAGATTCGCGAACATGCAGGATTGAGCTGT
 GCGAGGCATAGCATCCCCACATGGTAGACGAAGCGGATGTGGGTCTGATGGCTCTACTGCCGTCCAC
 CCGGGAGGATGGACGTGTGTTCCGGGATATCAGAGCTTGGCGTGACCAGCCACATGTCCATGCAGCC
 CTCCCGCTGGCTGTGTCGACATTTGGGTGTCACTGAAGGATTTACTGGTGCAGACTTTACAAAGGCA
 AAATTGTTTCTCAAGGAAGTGGACATGAAAGATCTAAAATCAGGGCTGCTTAAAAAAGCAAGCACAGA
 TTAATGACCTATCTACATCCTACACGGTTTAG

Restriction Sites: SgfI-MluI

ACCN: NM_001205273

Insert Size: 447 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.



[View online »](#)

RefSeq: NM_001205273.2
RefSeq Size: 974 bp
RefSeq ORF: 447 bp
Locus ID: 100507527
Cytogenetics: 19q12
MW: 16.1 kDa