

Product datasheet for SC329813

ATF3 (NM_001206488) Human Untagged Clone


Product data:

Product Type: Expression Plasmids
 Product Name: ATF3 (NM_001206488) Human Untagged Clone
 Tag: Tag Free
 Symbol: ATF3
 Vector: pCMV6-Entry (PS100001)
 Fully Sequenced ORF: >SC329813 representing NM_001206488.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGTCCTCTGCGCTGGAATCAGTCACTGTCAGCGACAGACCCCTCGGGGTGCCATCACAAAAGCCGAG
GTAGCCCTGAAGAAGATGAAAGGAAAAAGAGGCGACGAGAAAGAAATAAGATTGCAGCTGCAAAGTGC
CGAAACAAGAAGAAGGAGAAGACGGAGTGCCTGCAGAAAGAGTCGGAGAAGCTGGAAAGTGTGAATGCT
GAACTGAAGGCTCAGATTGAGGAGCTCAAGAACGAGAAGCAGCATTGATATACATGCTCAACCTTCAT
CGGCCACGTGTATTGTCCGGGCTCAGAATGGGAGGACTCCAGAAGATGAGAGAAACCTCTTTATCCAA
CAGATAAAGAAGGAACATTGCAGAGCTAA
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001206488

Insert Size: 375 bp



[View online >](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001206488.2
RefSeq Size:	2065 bp
RefSeq ORF:	375 bp
Locus ID:	467
UniProt ID:	P18847
Cytogenetics:	1q32.3
Protein Families:	Transcription Factors
MW:	14.3 kDa
Gene Summary:	<p>This gene encodes a member of the mammalian activation transcription factor/cAMP responsive element-binding (CREB) protein family of transcription factors. This gene is induced by a variety of signals, including many of those encountered by cancer cells, and is involved in the complex process of cellular stress response. Multiple transcript variants encoding different isoforms have been found for this gene. It is possible that alternative splicing of this gene may be physiologically important in the regulation of target genes. [provided by RefSeq, Apr 2011]</p> <p>Transcript Variant: This variant (8) differs in the 5' UTR, lacks a portion of the 5' coding region, and uses a downstream in-frame start codon, compared to variant 1. The encoded isoform (3) is shorter at the N-terminus, compared to isoform 1. Both variants 5 and 8 encode isoform 3.</p>