

Product datasheet for **SC329819**

TRIM49D1 (NM_001206627) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	TRIM49D1
Synonyms:	TRIM49D; TRIM49DIP; TRIM49D2; TRIM49D2P; TRIM49DP
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC329819 representing NM_001206627. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGAATTCTGGAATCTCGCAAGTCTTCCAGAGGGAACCTCACCTGCCCATCTGCCTGAACACTTTCATA
GACCCGGTCACCATAGACTGTGGGCACAGCTTTTGCAGGCCCTGTTTCTACCTCAACTGGCAAGACATC
CCAATTCTTACTCAGTCTTTGAATGCCTAAAGACAACACAGCAGAGAAACCTCAAACTAACATTTCGA
TTGAAGAAGATGGCTTCCCGTGCCAGAAAAGCCAGTCTCTGGCTATTCTGAGCTCTGAGGAGCAATG
TGTGGCACTCACAGGGAGACAAAGAAGATATTCTGTGAAGTGGACAGGAGCCTGCTCTGTTTGTCTGTC
TCCAGCTCTCTGGAGCACCGGTATCACAGACACTGTCCCGCTGAGTGGGCTGCTGAGGAACACCGGGAG
AAGCTTTTAAAGAAAAATGCAGTCTTTATGGGAAAAAGTTTGTGAAAATCAGAGAAACCTGAACGTGGAA
ACCACCAGAATCAGCCACTGGAAGGATTATGTGAATGTAAGGCTAGAAGCTATTAGAGCTGAGTATCAG
AAGATGCCTGCATTTTCATCATGAAGAAGAAAAACATAATTTGGAGATGCTGAAAAAGAAGGGGAAAGAT
ATTTTTCATCGACTTCATTTAAGTAAAGCCAAATGGCTCACAGGAGGGAGATTTTAAGAGGAACGTAT
GCGGAGCTGATGAAAAATGTCCATAAACCAGATGTGGAGCTACTTCAGGCTTTTGGAGACATATTACAC
AGGAGTGAGTCCGTGCTGCTGCACATGCCCAGCCTCTGAATCTAGAGCTCAGGGCAGGGCCCATCACT
GGACTGAGGGACAGGCTCAACCAATTCGAGTGGATATTACTCTGCCTCATAATGAAGCCAACAGTCAT
ATCTTCCGACGTGGAGATTTGAGAAGCATTGTATTGGATGTGACCGTCAAAATGCGCCCCATATCACT
GCAACACCTACAAGTTTCTTGATGGGGTGTCTCAGACTTTACCTCTGGCAAATATTACTGGGAGGTC
CATGTGGGGGACTCTTGAATTGGGCCTTTGGTGTCTGTAATAAGTATTGGAAGGGACGAATCAGAAT
GGCAATATACATGGAGAGGAGGAGTCTTTAGTCTTGGGTGTGTTAAGAATGACATTCACTGCAGCCTC
TTTACCACCTCCCCACTTACACTGCAGTATGTCCCAAGACCTACCAACCATGTAGGATTATTCCTGGAT
TGTGAAGCTAGAACTGTGAGCTTCGTTGATGTTAATCAAAGCTCCCCTATACACACCATCCCTAATTGC
TCCTTCTCACCTCCTCTCAGGCCTATCTTTGCTGTGTTACCTCTGA

```

Restriction Sites:	SgfI-MluI
ACCN:	NM_001206627
Insert Size:	1359 bp



OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001206627.1
RefSeq Size:	1474 bp
RefSeq ORF:	1359 bp
Locus ID:	399939
UniProt ID:	C9JIS8
Cytogenetics:	11q14.3
MW:	52.5 kDa