

Product datasheet for **SC329850**

CTNND1 (NM_001206884) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CTNND1 (NM_001206884) Human Untagged Clone
Tag:	Tag Free
Symbol:	CTNND1
Synonyms:	BCDS2; CAS; CTNND; p120; p120(CAS); p120(CTN); P120CAS; P120CTN
Vector:	pCMV6-Entry (PS100001)




[View online »](#)

Fully Sequenced ORF: >SC329850 representing NM_001206884.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGGCCAACGGCACACTCACCCGCCGCATCAGAACGGCCGGTTTGTGGGCGATGCTGACCTTGAAAGA
CAGAAATTTTCAGATTTGAAACTCAACGGACCCAGGATCACAGTCACCTTCTATATAGCACCATCCCC
AGGATGCAGGAGCCGGGCGAGATTGTGGAGACCTACACGGAGGAGGATCCTGAGGGAGCCATGTCTGA
GTCTCTGTGGAGACCTCAGATGATGGGACCACTCGGCCACAGAGACCACGGTCAAGAAAGTAGTGAAG
ACTGTGACAACACGGACAGTACAGCCAGTCGCTATGGGACCAGACGGGTTGCCTGTGGATGCTTCATCA
GTTTCTAACAACTATATCCAGACTTTGGGTCGTGATTTCCGCAAGAATGGCAATGGGGGACCTGGTCCC
TATGTGGGGCAAGCTGGCACTGCTACCCTTCTAGGAACTTCCACTACCCTCCTGATGGTTATAGTCGC
CACTATGAAGATGGTTATCCAGGTGGCAGTGATAACTATGGCAGTCTGTCCCGGGTGACCCGATTGAG
GAGCGGTATAGGCCCAGCATGGAAGGCTACCGGGCACCTAGTAGACAGGATGTGTATGGGCCCAACCC
CAGGTTCCGGTAGGTGGGAGCAGCGTGGATCTGCATCGTTTTATCCAGAGCCTTATGGGCTAGAGGAT
GACCAGCGTAGTATGGGCTATGATGACCTGGATTATGGTATGATGTCTGATTATGGCACTGCCCGTCGG
ACTGGGACACCCTCTGACCCTCGTCGGCGCCTCAGGAGCTATGAAGACATGATTGGTGGAGAGGTGCCA
TCGGATCAATACTACTGGGCTCCTTTGGCCAGCATGAGCGAGGAAGTTTAGCAAGCTTGGATAGCCTG
CGCAAAGGAGGGCCTCCACCTCCTAATTGGAGACAGCCAGAGCTGCCAGAGGTGATCGCCATGCTTGG
TTCCGCTTGGATGCTGTCAAGTCCAATGCAGCTGCATACCTGCAACACTTATGCTACCGCAATGACAAG
GTGAAGACTGACGTGCGGAAGCTCAAGGGCATCCCAGTACTGGTGGGATTGTTAGACCATCCCAAAAAG
GAAGTGACCTTGGAGCCTGTGGAGCTCTCAAGAAATATCTTTTTGGACGTGACCAGGATAACAAGATT
GCCATAAAAACTGTGATGGTGTGCCTGCCCTTGTGCGATTGCTTCGAAAGGCTCGTGATATGGACCTT
ACTGAAGTTATTACCGGAACCTGTGGAATCTTTCATCCCAGTCAATCAAAAATGGAGATTGTGGAC
CATGCACTGCATGCCTTGACAGATGAAGTGAATTCCTCATTCTGGTTGGGAGCGGGAACCTAATGAA
GACTGTAAAGCCACGCCATATTGAGTGGGAATCGGTGCTACCAACACAGCTGGCTGCCTTAGGAATGTA
AGCTCAGAGAGGAGTGAAGCTCGCCGAAACTCGGGAATGTGATGGTTTAGTTGATGCCCTCATTTTC
ATTGTTCAAGGCTGAGATTGGGCAGAAGGATTAGAGCAAGCTTGTAGAGAAGTGTGTTTGCCTTCTT
CGGAACTTATCATATCAAGTTCACCGGGAGATCCACAGGCAGAGCGTTACCAAGAGGCAGCTCCCAAT
GTTGCCAACAATACTGGGCCACATGCTGCCAGTTGCTTTGGGGCCAAGAAGGGCAAAGATGAGTGGTTC
TCCAGAGGGAAAAACCTATAGAGGATCCAGCAAACGATACAGTGGATTTCCCTAAAAGAACGAGTCCA
GCTCGAGGCTATGAGCTTTATTTAGCCAGAGGTGGTTCGGATATACATCTCACTTCTTAAGGAGAGC
AAGACTCCTGCCATCCTAGAAGCCTCAGCTGGAGCTATCCAGAAGTGTGTGCTGGGCGCTGGACGTAT
GGTCGATACATCCGCTCTGCTCTGCGTCAAGAGAAGGCTTTTCTGCCATAGCTGACCTCCTGACTAAT
GAACATGAACGGGTGGTAAAGCTGCATCTGGAGCACTGAGAAACCTGGCTGTGGATGCTCGCAACAAA
GAATTAATTGGTAAACATGCTATTCTAACTTGGTAAAGAATCTGCCAGGAGGACAGCAGAAGTCTCTT
TGGAATTTCTCTGAGGACACTGTCACTCTATTTTGAACACTATCAACGAGGTTATCGCTGAGAAGTTG
GAGGCTGCCAAAAAGCTTCGAGAGACACAGGATTTGAGAAGCTGGTGTGATCAACAAATCAGGGAAC
CGCTCAGAAAAAGAAGTTCGAGCAGCAGCACTTGATTACAGACAATCTGGGGATATAAGGAACTGCGG
AAGCCACTGAAAAAGAAGGATGGAAGAAATCAGACTTTCAGGTGAATCTAAACAATGCTTCCCGAAGC
CAGAGCAGTCATTATATGATGATGATGATGATGATGATGATGATGATGATGATGATGATGATGATGAT
CCTGATCGGGAAGAAATTCAGATGAGCAATATGGGATCAAACACAAAATCACTAGATAAACAATTTCC
ACACCAAATGAGAGAGGAGACCACAATAGAACACTGGATCGATCGGGGGTCTAGGCGACATGGAGCCA
TTGAAGGGAACAACACCTTGTATGCAGAAGATTAG
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_001206884

Insert Size: 2658 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001206884.1](#)

RefSeq Size: 5975 bp

RefSeq ORF: 2658 bp

Locus ID: 1500

UniProt ID: [O60716](#)

Cytogenetics: 11q12.1

Protein Families: Druggable Genome

Protein Pathways: Adherens junction, Leukocyte transendothelial migration

MW: 98.9 kDa

Gene Summary: This gene encodes a member of the Armadillo protein family, which function in adhesion between cells and signal transduction. Multiple translation initiation codons and alternative splicing result in many different isoforms being translated. Not all of the full-length natures of the described transcript variants have been determined. Read-through transcription also exists between this gene and the neighboring upstream thioredoxin-related transmembrane protein 2 (TMX2) gene. [provided by RefSeq, Dec 2010]
Transcript Variant: This variant (15) differs in the 5' UTR and has multiple coding region differences, compared to variant 1. The encoded isoform (2AC) is shorter than isoform 1ABC.
Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.