

## Product datasheet for **SC329856**

### CTNND1 (NM\_001206890) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	CTNND1 (NM_001206890) Human Untagged Clone
Tag:	Tag Free
Symbol:	CTNND1
Synonyms:	BCDS2; CAS; CTNND; p120; p120(CAS); p120(CTN); P120CAS; P120CTN
Vector:	pCMV6-Entry (PS100001)




[View online »](#)

Fully Sequenced ORF: >SC329856 representing NM\_001206890.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGCAGGAGCCGGGGCAGATTGTGGAGACCTACACGGAGGAGGATCCTGAGGGAGCCATGTCTGTAGTC
TCTGTGGAGACCTCAGATGATGGGACCACTCGGCCACAGAGACCACGGTCAAGAAAGTAGTGAAGACT
GTGACAACACGGACAGTACAGCCAGTCGCTATGGGACCAGACGGGTTGCCTGTGGATGCTTCATCAGTT
TCTAACAACTATATCCAGACTTTGGGTCGTGATTTCCGCAAGAATGGCAATGGGGGACCTGGTCCCTAT
GTGGGGCAAGCTGGCACTGCTACCCTTCTAGGAACTTCCAATACTACCCTCCTGATGGTTATAGTCGCCAC
TATGAAGATGGTTATCCAGGTGGCAGTGATAACTATGGCAGTCTGTCCCGGGTGACCCGCATTGAGGAG
CGGTATAGGCCACAGATGGAAGGCTACCGGGCACCTAGTAGACAGGATGTGTATGGGCCCAACCCAG
GTTCCGGTAGGTGGGAGCAGCGTGGATCTGCATCGCTTTCATCCAGAGCCTTATGGGCTAGAGGATGAC
CAGCGTAGTATGGGCTATGATGACCTGGATTATGGTATGATGTCTGATTATGGCACTGCCCGTCGGACT
GGGACACCCTCTGACCCTCGTCGGCGCCTCAGGAGCTATGAAGACATGATTGGTGAGGAGGTGCCATCG
GATCAATACTACTGGGCTCCTTTGGCCAGCATGAGCGAGGAAGTTTAGCAAGCTTGGATAGCCTGCCG
AAAGGAGGGCCTCCACCTCCTAATTGGAGACAGCCAGAGCTGCCAGAGGTGATCGCCATGCTTGGATTC
CGCTTGGATGCTGTCAAGTCCAATGCAGCTGCATACCTGCAACACTTATGCTACCCAATGACAAGGTG
AAGACTGACGTGCGGAAGCTCAAGGGCATCCCAGTACTGGTGGGATTGTTAGACCATCCCAAAAAGGAA
GTGCACCTTGGAGCCTGTGGAGCTCTCAAGAATATCTCTTTTGGACGTGACCAGGATAACAAGATTGCC
ATAAAAACTGTGATGGTGTGCCTGCCCTTGTGCGATTGCTTCGAAAGGCTCGTGATATGGACCTTACT
GAAGTTATTACCGGAACCTGTGGAATCTTTCATCCCATGACTCAATCAAAATGGAGATTGTGGACCAT
GCACTGCATGCCTTGACAGATGAAGTGATCATTCTCATTCTGGTTGGGAGCGGGAACTAATGAAGAC
TGTAAGCCACGCCATATTGAGTGGGAATCGGTGCTCACCAACACAGCTGGCTGCCTTAGGAATGTAAGC
TCAGAGAGGAGTGAAGCTCGCCGAACTTCGGGAATGTGATGGTTTAGTTGATGCCCTCATTTTCATT
GTTCAGGCTGAGATTGGGCAGAAGGATTACAGACAGCAAGCTTGTAGAGAAGTGTGTTTGCCTTCTCGG
AACTTATCATATCAAGTTCACCGGGAGATCCACAGGCAGAGCGTTACCAAGAGGAGCTCCCAATGTT
GCCAACATACTGGGCCACATGCTGCCAGTTGCTTTGGGGCAAGAAGGGCAAAGGGAAAAAACCTATA
GAGGATCCAGCAAACGATACAGTGGATTCCCTAAAAGAACGAGTCCAGCTCGAGGCTATGAGCTCTTA
TTTCAGCCAGAGGTGGTTCGGATATACATCTCACTTCTTAAGGAGAGCAAGACTCCTGCCATCCTAGAA
GCCTCAGCTGGAGCTATCCAGAAGTGTGTGCTGGGCGCTGGACGTATGGTCGATACATCCGCTCTGCT
CTGCGTCAAGAGAAGGCTCTTTCTGCCATAGCTGACCTCCTGACTAATGAACATGAACGGGTGGTAAA
GCTGCATCTGGAGCACTGAGAACTGGCTGTGGATGCTCGCAACAAGAATTAATTGGTAAACATGCT
ATTCCTAACTTGGTAAAGAATCTGCCAGGAGGACAGCAGAAGCTCCTTGGAAATTTCTCTGAGGACT
GTCATCTCTATTTTGAACACTATCAACGAGGTTATCGCTGAGAAGTGGAGGCTGCCAAAAAGCTTCGA
GAGACACAGGGTATTGAGAAGCTGGTGTGATCAACAAATCAGGGAACCGCTCAGAAAAAGAAGTTCGA
GCAGCAGCACTTGTATTACAGACAATCTGGGGATATAAGGAAGTCCGGAAGCCACTGGAAAAAGAAGGA
TGAAGAAATCAGACTTTCAGGTGAATCTAAACAATGCTTCCCGAAGCCAGAGCAGTCATTTCATATGAT
GATAGTACTCTCCCTCTCATTGACCGGAACCAAAAATCAGATAAGAAACCTGATCGGGAAGAAATTCAG
ATGAGCAATATGGGATCAACACAAAATCACTAGATAACAACATTTCCACACCAATGAGAGAGGAGAC
CACAAATAGAACACTGGATCGATCGGGGATCTAGGCGACATGGAGCCATTGAAGGGAACAACACCCCTTG
ATGCAGAAGATTAG
```

Restriction Sites: SgfI-MluI

## Plasmid Map:



ACCN: NM\_001206890

Insert Size: 2499 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001206890.1](#)

RefSeq Size: 5799 bp

RefSeq ORF: 2499 bp

Locus ID: 1500

UniProt ID: [O60716](#)

Cytogenetics: 11q12.1

Protein Families: Druggable Genome

Protein Pathways: Adherens junction, Leukocyte transendothelial migration

**MW:** 92.7 kDa

**Gene Summary:** This gene encodes a member of the Armadillo protein family, which function in adhesion between cells and signal transduction. Multiple translation initiation codons and alternative splicing result in many different isoforms being translated. Not all of the full-length natures of the described transcript variants have been determined. Read-through transcription also exists between this gene and the neighboring upstream thioredoxin-related transmembrane protein 2 (TMX2) gene. [provided by RefSeq, Dec 2010]  
Transcript Variant: This variant (21) differs in the 5' UTR and has multiple coding region differences, compared to variant 1. The encoded isoform (3A) is shorter than isoform 1ABC.  
Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.