

Product datasheet for **SC329886**

ETV7 (NM_001207035) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: ETV7 (NM_001207035) Human Untagged Clone
Tag: Tag Free
Symbol: ETV7
Synonyms: TEL-2; TEL2; TELB
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC329886 representing NM_001207035.
Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
ATGCAGGAGGGAGAATTGGCTATTTCTCCTATAAGCCCTGTGGCAGCCATGCCTCCCCTAGGCACCCAC
GTGCAAGCCAGATGTGAAGCTCAAATTAACCTGCTGGGTGAAGGGGGATCTGCAAGCTGCCAGGAAGA
CTCCGCATCCAGCCCGCACTGTGGAGCAGGGAGGACGTGCTGCACTGGCTGCGCTGGCAGAGCAGGAG
TACTCTCTGCCATGCACCGGGAGCACGGGTTTCGAGATGAACGGACGCGCCCTCTGCATCCTCACCAAG
GACGACTTCCGGCACCGTGGCCAGCTCAGGTGACGTCCTGTATGAGCTGCTCCAGTACATCAAGACC
CAGCGGGAGCCCTGGTGTGTGGCCCTTTTTTGGAGGGATCTTCAGGCTGAAGACGCCACCCAGCAC
TCTCCAGTCCCCCGGAAGAGGTGACTGGCCCTCTCAGATGGACACCCGAAGGGGCCACCTGCTGCAG
CCACCAGACCCAGGGCTTACCAGCAACTTCGGCCACCTGGATGACCCTGGCTGGCAAGGTGGACCCCT
GGCAAGGAGGAGTCCCTCAACTTATGCTACTGTGCAGAGCTCGGCTGCAGGACCCAGGGGGTCTGTTCC
TTCCCCGCGATGCCCGAGGCCCCATTGACGGCAGGATCGCTGACTGCCGCTGCTGTGGGATTACGTG
TATCAGCTGCTCCTTGATACCCGATATGAGCCCTACATCAAGTGGGAAGACAAGGACGCCAAGATCTTC
CGAGTTGTGGATCCAAATGGGCTCGCCAGACTCTGGGAAATCACAAGAACCGGGTGAACATGACCTAC
GAGAAGATGTCTCGTGCCTGCGCCACTATTATAAGCTTAATATCATTAAAGAAGAACCGGGGCAGAAA
CTCCTGTTCAGAAATGGACTTCAGCTGATCTTCATATTCATATGGAGTTTCCAGTGA
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:



ACCN: NM_001207035

Insert Size: 954 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq:	NM_001207035.1
RefSeq Size:	1634 bp
RefSeq ORF:	954 bp
Locus ID:	51513
UniProt ID:	Q9Y603
Cytogenetics:	6p21.31
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Dorso-ventral axis formation
MW:	36.2 kDa
Gene Summary:	<p>The protein encoded by this gene belongs to the ETS family of transcription factors, which is a large group of evolutionarily conserved transcriptional regulators that play an important role in a variety of cellular processes throughout development and differentiation, and are involved in oncogenesis as well. This protein is predominantly expressed in hematopoietic tissues. Several alternatively spliced transcript variants encoding different isoforms have been described for this gene (PMID:11108721).[provided by RefSeq, May 2011]</p> <p>Transcript Variant: This variant (2) contains an alternate 3' terminal exon compared to variant 1. This results in a frame-shift and a shorter isoform (2, also known as Tel-2e) with a distinct C-terminus compared to isoform 1.</p>