

Product datasheet for **SC329901**

CXADR (NM_001207063) Human Untagged Clone

Product data:


Product Type: Expression Plasmids
Product Name: CXADR (NM_001207063) Human Untagged Clone
Tag: Tag Free
Symbol: CXADR
Synonyms: CAR; CAR4/6; HCAR
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC329901 representing NM_001207063.
Blue=Insert sequence **Red**=Cloning site **Green**=Tag(s)

```
ATGGCGCTCCTGCTGTGCTTCGTGCTCCTGTGCGGAGTAGTGGATTCGCCAGAAGTTTGAGTATCACT
ACTCCTGAAGAGATGATTGAAAAAGCCAAAGGGGAACTGCCTATCTGCCATGCAAATTTACGCTTAGT
CCCGAAGACCAGGGACCCTGGACATCGAGTGGCTGATATCACCAGCTGATAATCAGAAGGTGGATCAA
GTGATTATTTTATTCTGGAGACAAAATTTATGATGACTACTATCCAGATCTGAAAGGCCGAGTACAT
TTTACGAGTAATGATCTCAAATCTGGTATGCATCAATAAATGTAAACGAATTTACAAGTGCAGATATT
GGCACATATCAGTGCAAAGTGAAAAAGCTCCTGGTGTGCAAAATAAGAAGATTCATCTGGTAGTTCTT
GTTAAGCCTTCAGGTGCGAGATGTTACGTTGATGGATCTGAAGAAATTTGGAAGTGACTTTAAGATAAAA
TGTGAACCAAAAGAAGGTTCACTTCCATTACAGTATGAGTGGCAAAAATTGTCTGACTCACAGAAAATG
CCCACTTCATGGTTAGCAGGGAAGATGTGCCACCTCCAAGAGCCGTACGTCCACTGCCAGAAGCTACA
TCGGCAGTAATCATTATCCCTGGGGTCCATGTCTCCTTCCAACATGGAAGGATATTCCAAGACTCAGT
ATAACCAAGTACCAAGTGAAGACTTTGAACGCACTCCTCAGAGTCCGACTCTCCACCTGCTAAGGTAG
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:


ACCN: NM_001207063

Insert Size: 759 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001207063.1</u>
RefSeq Size:	5482 bp
RefSeq ORF:	759 bp
Locus ID:	1525
UniProt ID:	<u>P78310</u>
Cytogenetics:	21q21.1
Protein Families:	Druggable Genome, Stem cell - Pluripotency, Transmembrane
Protein Pathways:	Viral myocarditis
MW:	28.2 kDa
Gene Summary:	<p>The protein encoded by this gene is a type I membrane receptor for group B coxsackieviruses and subgroup C adenoviruses. Several transcript variants encoding different isoforms have been found for this gene. Pseudogenes of this gene are found on chromosomes 15, 18, and 21. [provided by RefSeq, May 2011]</p> <p>Transcript Variant: This variant (2, alternatively referred to as beta or CAR4/7) lacks two consecutive alternate coding exons compared to variant 1, that causes a frameshift. The resulting isoform (2) has a shorter and distinct C-terminus compared to isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>