

Product datasheet for SC329920

USP17L24 (NM_001242327) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	USP17L24 (NM_001242327) Human Untagged Clone
Tag:	Tag Free
Symbol:	USP17L24
Synonyms:	USP17L25; USP17L26; USP17L27; USP17L28; USP17L29; USP17L30
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC329920 representing NM_001242327. Blue=Insert sequence Red=Cloning site Green=Tag(s)

ATGGAGGACGACTCACTCTACTTGAGAGGTGAGTGGCAGTTCAACCACTTTTCAAACACTCACATCTTCT
 CGGCCCGATGCAGCTTTTCTGAAATCCAGCGGACTTCTCTCCCTGAGAAGTCACCACTCTCATGTGAG
 ACCCGTGTGACCTCTGTGATGATTGGCTCCTGTGGCAAGACAGCTTGCTCCAGGGAGAAGCTTCTCT
 CTGAGTAGCAGGAGACCTGCTGCGGTGGGGCTGGGCTCCAGAATATGGGAAATACCTGCTACGTGAAC
 GCTTCTTGCAGTGCCTGACATACACACCGCCCTTGCCAACTACATGCTGTCCCGGAGCACTCTCAA
 ACGTGTCACTCGTCAAGGGCTGCATGCTGTACTATGCAAGCTCACATCACACGGGCCCTCCACAAT
 CCTGGCCACGTATCCAGCCCTCACAGGCATTGGCTGCTGGCTTCCATAGAGGCAAGCAGGAAGATGCC
 CATGAATTTCTCATGTTCACTGTGGATGCCATGAAAAAGGCATGCCTTCCCGGGCACAAGCAGGTAGAT
 CATCACTCTAAGGACACCACTCATCCACCAATATTTGGAGGCTACTGGAGATCTCAAATCAAGTGT
 CTCCACTGCCACGGCATTTCAGACACTTTTGACCCTTACCTGGACATCGCCCTGGATATCCAGGCAGCT
 CAGAGTGTCCAGCAAGCTTTGGAACAGTTGGTGAAGCCCGAAGAAGCTCAATGGAGAGAATGCCTATCAT
 TGTGGTGTGTGCTCCAGAGGGCGCCGGCTCCAAGACGTTAACTTTACACACCTCTGCCAAGGTCTCTC
 ATCCTTGATTGAAGAGATTCTCCGATGTCACAGGCAACAAGATTGCCAAGAATGTGCAATATCCTGAG
 TGCCTTGACATGCAGCCATACATGTCTCAGCCGAACACAGGACCTCTCGTCTATGCTCTATGCTGTG
 CTGGTCCACGCTGGGTGGAGTTGTCACAACGGACATTACTTCTCTTATGTCAAAGCTCAAGAAGGCCAG
 TGGTATAAAATGGATGATGCCGAGGTCAACGCTCTAGCATCACTTCTGTCTGAGTCAACAGGCCTAC
 GTCCTCTTTTACATCCAGAAGAGTGAATGGGAAAGACACAGTGAAGTGTGTCAAGAGGCAGGGAACCA
 AGAGCCCTTGGCGCAGAAGACACCGACAGGCGAGCAACGCAAGGAGAGCTCAAGAGAGACCACCCCTGC
 CTCCAGGCCCCGAGTTGGACGAGCACTTGGTGGAAAGAGCCACTCAGGAAAGCACCTTAGACCACTGG
 AAATTCCTTCAAGAGCAAAACAAAACGAAGCCTGAGTTCAACGTCAGAAAAGTCGAAGGTACCCTGCCT
 CCCGACGTAATTTGATTCAATCAAAATACAAGTGTGGGATGAAGAACCATCATCCTGAACAGCAA
 AGCTCCCTGCTAAACCTCTCTTCGTCGACCCCGACACATCAGGAGTCCATGAACACTGGCACACTCGCT
 TCCTGCGAGGGAGGCCAGGAGATCAAAGGGAAGAACAACACAGCAAGAGGGCTCTGCTTGTGTGC
 CAGTGA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001242327


[View online »](#)

Insert Size:	1593 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001242327.1</u>
RefSeq Size:	1593 bp
RefSeq ORF:	1593 bp
Locus ID:	728369
UniProt ID:	<u>Q0WX57</u>
Cytogenetics:	4p16.1
MW:	59.7 kDa