

Product datasheet for SC329921

USP17L26 (NM_001242328) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: USP17L26 (NM_001242328) Human Untagged Clone
Tag: Tag Free
Symbol: USP17L26
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC329921 representing NM_001242328.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGAGGACGACTCACTCTACTTGAGAGGTGAGTGGCAGTTCAACCACCTTTCAAACACTCACATCTTCT
CGGCCCGATGCAGCTTTTGGCTGAAATCCAGCGGACTTCTCTCCCTGAGAAGTCACCACCTCATGTGAG
ACCCGTGTCGACCTCTGTGATGATTTGGCTCCTGTGGCAAGACAGCTTGTCCCAGGGAGAAGCTTCTCT
CTGAGTAGCAGGAGACCTGCTGCGGTGGGGCTGGGCTCCAGAATATGGGAAATACCTGTACGTGAAC
GCTTCTTGCAGTGCCTGACATACACACCGCCCTTCCCAACTACATGCTGTCCCGGGAGCACTCAA
ACGTGTCATCGTCACAAGGGCTGCATGCTGTACTATGCAAGTCACATCACACGGCCCTCCACAAT
CCTGGCCACGTATCCAGCCCTCACAGGCAATTGGCTGCTGGCTTCCATAGAGGCAAGCAGGAAGATGCC
CATGAATTTCTCATGTTCACTGTGGATGCCATGAAAAAGGCATGCCTTCCCGGGCACAAGCAGGTAGAT
CATCACTCTAAGGACACCACCCTCATCCACCAAATATTTGGAGGCTACTGGAGATCTCAAATCAAGTGT
CTCCACTGCCACGGCATTTCAGACACTTTTGACCCTTACCTGGACATCGCCCTGGATATCCAGGCAGCT
CAGAGTGTCCAGCAAGCTTTGGAACAGTTGGTGAAGCCGGAAGAAGCACTCAATGGAGAGAATGCCTATCAT
TGTGGTGTGTCTCCAGAGGGCGCCGGCCTCAAAGACGTTAACTTTACACACCTCTGCCAAGGTCTC
ATCCTTGTATTGAAGAGATTCTCCGATGTCACAGGCAACAAGATTGCCAAGAATGTGCAATATCCTGAG
TGCCTTGACATGCAGCCATACATGTCTCAGCCGAACACAGGACCTCTCGTCTATGCTCTATGCTGTG
CTGGTCCACGCTGGGTGGAGTTGTACAAACGGACATTACTTCTCTATGTCAAAGCTCAAGAAGGCCAG
TGGTATAAAATGGATGATGCCGAGGTACCCGCCTTAGCATCACTTCTGTCTGAGTCAACAGGCCCTAC
GTCCTCTTTTACATCCAGAAGAGTGAATGGGAAAGACACAGTGAGAGTGTGTCAAGAGGCAGGGAACCA
AGAGCCCTTGGCGCAGAAGACACCGACAGGCGAGCAACGCAAGGAGAGCTCAAGAGAGACCACCCCTGC
CTCCAGGCCCCCGAGTTGGACGAGCACTTGGTGGAAAGAGCCACTCAGGAAAGCACCTTAGACCCTGG
AAATTCCTTCAAGAGCAAAACAAAACGAAGCCTGAGTTCAACGTCAGAAAAGTCGAAGGTACCCTGCCT
CCCAGCTACTTGTGATTCATCAATCAAATACAAGTGTGGATGAAGAACCATCATCCTGAACAGCAA
AGCTCCCTGCTAAACCTCTCTTCGTCGACCCCGACACATCAGGAGTCCATGAACACTGGCACACTCGCT
TCCCTGCCGAGGGAGGCCAGGAGATCCAAGGGAAGAACAACACAGCAAGAGGGCTCTGCTTGTGTGC
CAGTGA
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001242328
Insert Size: 1593 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001242328.1</u>
RefSeq Size:	1593 bp
RefSeq ORF:	1593 bp
Locus ID:	728379
UniProt ID:	<u>Q0WX57</u>
Cytogenetics:	4p16.1
MW:	59.7 kDa