

Product datasheet for **SC329967**

FAM156A (NM_001242493) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FAM156A (NM_001242493) Human Untagged Clone
Tag:	Tag Free
Symbol:	FAM156A
Synonyms:	PRO0659; TMEM29
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC329967 representing NM_001242493. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
ATGGATCCA CTCCAGAAACGGAATCCAGCATCGCCTTCCAATCTTCCCGATGACAGCTGCAGAGACT
TCCCAGGAAGGTCCAGCGCCCTCTCAGCCTTCGTACTCAGAACAGCCGATGATGGGCCTCAGTAACCTG
AGCCCCGGTCCTGGCCCCAGCCAGGCCGTGCCTCTCCAGAGGGGCTGCTCCGCCAGCGGTACAGAGAG
GAGAAGACCCTGGAAGAGCGGCCGTGGGAGAGGCTGGAGTTCCTTCAGAGGAAGAAAGCATTCTGCCG
CATGTGAGGAGGAGACACCGCATCACATGGCCCCCTATGCTGTTGGGAGGGAAGCCAGAATCTCCCA
TTAGGTGACAGAAGTCAGAATCGATTCCGATGTGAATGTCGATACTGCCAGAGCCACAGGCCGAATCTT
TCTGGGATCCCTGGGAGAGTAACAGGGCCCCACATCCCTCCTCTGGGAGACGCTGGTGCAGGGCCTC
AGTGGCTTGACTCTCAGCCTAGGCACCAACCAGCCGGGCTCTGCCTGAAGCGGCACTCCAGCCACAG
GAGACAGAGGAGAAGCGCCAGCGAGAGAGGCAGCAGGAGAGCAAATAATGTTTCAGAGGCTGCTCAAG
CAGTGGTTAGAGGAAAACTGA
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:



ACCN: NM_001242493

Insert Size: 642 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001242493.1](#)

RefSeq Size: 1839 bp

RefSeq ORF: 642 bp

Locus ID: 29057

UniProt ID: [Q8NDB6](#)

Cytogenetics: Xp11.22

MW: 24.4 kDa