

Product datasheet for **SC329973**

MIF4GD (NM_001242500) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	MIF4GD
Synonyms:	AD023; MIFD; SLIP1
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC329973 representing NM_001242500. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGGGGAGCCAGTAGAGAGGAGTATAAAATCCAGTCCTTTGATGCAGAGACCCAGCAGCTGCTGAAG
ACAGCACTCAAAGTTGCCTGCTTTGAGACAGAAGATGGAGAATATTCGTTTGCCAGAGGAGCTATAGT
AATTGCTCTCGCTGATGCCTTCCCGCTGTAACACGCAGTACAGAGATCCGGGTGCTGTGGACTTGGAG
AAAGTGCCAATGTGATTGTGGACCATTCTCTGCAGGACTGTGTTCAGCAAGGAAGCAGGACGCATG
TGCTACGCCATCATTACGGCAGAGAGTAAACAAGCAGGCCAGAGTGTCTTCGACGTGGACTCCTCAAC
CGGCTGCAGCAGGAGTACCAGGCTCGGGAGCAGCTGCGAGCACGCTCCCTGCAGGGCTGGGTCTGCTAT
GTCACCTTTATCTGCAACATCTTTGACTACCTGAGGGTGAACAACATGCCCATGATGGCCCTGGTGAAC
CCTGTCTATGACTGCCTCTTCCGGCTGGCCAGCCAGACAGTTTGAGCAAGGAGGAGGAGGTGACTGT
TTGGTGTCTGAGCTGCACCGGTTGGGGAGCAGCTGGAGAAAATGAATGGGCAGCGCATGGATGAGCTC
TTTGTGCTGATCCGGGATGGCTTCTGCTCCCAACTGGCCTCAGCTCCCTGGCCAGCTGCTGCTGCTG
GAGATCATTGAGTTCGGGGCGGCCGCTGGAAGACAACGCCAGCTGCCACAAGTATTACTACAGCGAA
GTCTCCGACTAG

```

Restriction Sites:	SgfI-MluI
ACCN:	NM_001242500
Insert Size:	771 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001242500.1</u>
RefSeq Size:	1364 bp
RefSeq ORF:	771 bp
Locus ID:	57409
UniProt ID:	<u>A9UHW6</u>
Cytogenetics:	17q25.1
MW:	29.4 kDa
Gene Summary:	This gene encodes a protein which interacts with the N-terminus of the stem-loop binding protein (SLBP) and the 3' end of histone mRNA. This interaction facilitates the activation of histone mRNA translation. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jun 2011] Transcript Variant: This variant (3), as well as variant 2, encodes isoform (2).