

Product datasheet for **SC329977**

GDA (NM_001242505) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: GDA (NM_001242505) Human Untagged Clone
Tag: Tag Free
Symbol: GDA
Synonyms: CYPIN; GAH; GUANASE; NEDASIN
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC329977 representing NM_001242505.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

ATGTGTGCCGCTCAGATGCCGCCCTGGCGCACATCTTCCGAGGGACGTTTCGTCCACTCCACCTGGACC
TGCCCCATGGAGGTGCTGCGGGATCACCTCCTCGGCGTGAGCGACAGCGGCAAAATAGTGTTTTTAGAA
GAAGCATCTCAACAGGAAAACTGGCCAAAGAATGGTGCCTCAAGCCGTGTGAAATAAGAGAAGTGGC
CACCATGAGTTCTTCATGCCTGGGCTGGTTGATACACACATCCATGCCTCTCAGTATTCCTTTGCTGGA
AGTAGCATAGACCTGCCACTCTGGAGTGGCTGACCAAGTACACATTTCCCTGCAGAACACAGATTCCAG
AACATCGACTTTGCAGAAGAAGTATATACCAGAGTTGTCAGGAGAACACTAAAGAATGGAACAACCACA
GCTTGTTACTTTGCAACAATTCACACTGACTCATCTCTGCTCCTTGCCGACATTACAGATAAATTTGGA
CAGCGGGCATTGTGGGCAAAGTTTGCATGGATTTGAATGACACTTTTCCAGAATACAAGGAGACCACT
GAGGAATCGATCAAGGAAACTGAGAGATTTGTGTGAGAAATGCTCCAAAAGAACTATTCTAGAGTGAAG
CCCATAGTGACACCACGTTTTTCCCTCCTCTGCTGAGACTTTGATGGGTGAACTGGGCAACATTGCT
AAAACCCGTGATTTGCACATTGAGAGCCATATAAGTAAAATCGTGATGAAGTTGAAGCTGTGAAAAAC
TTATACCCAGTTATAAAAACTACACATCTGTGTATGATAAAAAACAATCTTTGACAAATAAGACAGTG
ATGGCACACGGCTGCTACCTCTCTGCAGAAGAACTGAACGTATTCCATGAACGAGGAGCATCCATCGCA
CACTGTCCCAATTCTAATTTATCGCTCAGCAGTGGATTTCTAAATGTGCTAGAAGTCTGAAACATGAA
GTCAAGATAGGGCTGGGTACAGACGTGGCTGGTGGCTATTTCATATTCCATGCTTGATGCAATCAGAAGA
GCAGTGATGGTTTCCAATATCCTTTTAATTAATAAGGTAATGAGAAAAGCCTCACCCCTCAAAGAAGTC
TTCAGACTAGCTACTCTTGGAGGAAGCCAAGCCCTGGGGCTGGATGGTGAGATTGAAAACCTTTGAAGTG
GGCAAGGAATTTGATGCCATCCTGATCAACCCAAAGCATCCGACTCTCCATTGACCTGTTTTATGGG
GACTTTTTTGGTGATATTTCTGAGGCTGTTATCCAGAAGTTCTCTATCTAGGAGATGATCGAAATATT
GAAGAGGTTTATGTGGGCGGAAAGCAGGTGGTTCCGTTTTCCAGCTCAGTGAAGAAACAATTCACTTA
CCAGCCTCCTCACCCATCCTCCACCATTTCCTTAA
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001242505

Insert Size: 1416 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001242505.2](#)

RefSeq Size: 2096 bp

RefSeq ORF: 1416 bp

Locus ID: 9615

UniProt ID: [Q9Y2T3](#)

Cytogenetics:	9q21.13
Protein Pathways:	Metabolic pathways, Purine metabolism
MW:	52.8 kDa
Gene Summary:	<p>This gene encodes an enzyme responsible for the hydrolytic deamination of guanine. Studies in rat ortholog suggest this gene plays a role in microtubule assembly. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2011]</p> <p>Transcript Variant: This variant (1) encodes the longest isoform (a). Variants 1 and 6 encode the same isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>