

Product datasheet for **SC329981**

MFSD11 (NM_001242533) Human Untagged Clone

Product data:


Product Type: Expression Plasmids
Product Name: MFSD11 (NM_001242533) Human Untagged Clone
Tag: Tag Free
Symbol: MFSD11
Synonyms: ET
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC329981 representing NM_001242533.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGTCCCCGGAATCTAAAAAGCTTTTCAACATCATTATTTTAGGAGTTGCCTTTATGTTTATGTTCACT
GCCTTTCAAAGTGTGGAAATGTGGCGCAAAGTGCATCAGGAGCTTAAATAGGACAGATTTTCACGGC
AGTGGATATACCAGCATGGCTATTATCTATGGAGTGTCTCTGCTTCAAATTTGATTACACCGTCAGTG
GTTGCCATTGTAGGACCTCAACTCTCTATGTTTCCAGTGGTTTATTTTACAGCATGTACATTGCCGTT
TTCATCCAGCCTTCCCCTGGTCTTCTACACAGCCTCTGTTTTCATTGGAATTGCTGCTGCTGTGCTT
TGGACAGCACAAGGAAAGTGCCTGACAATCAATTCGGATGAGCACAGCATTGGGAGAAACAGTGGGATT
TTCTGGGCACTTCTGCAGTCTAGCTTGTCTTTGGAAATCTCTACATATATTTTGCCTGGCAAGGGAAA
ACTCAGATATCAGAGAGTGACCGAAGAAGAGTGTATTGCCCTAACGGTGATTAGCCTTGTGGGGACA
GTTCTATTCTTTCTATTTCGAAACCAGATTCTGAAAATGTCCTAGGAGAAGATGAGTCTTCTGATGAC
CAGGACATGGAAGTCAACGAGTCTGCCCAGAACAATCTGACAAAGGCAGTAGATGCTTTTAAAAAGTCT
TTTAAGTTATGTGTCACCAAGGAGATGCTCCTTCTTAGTATTACAAGTCTTATACAGGTCTGGAATTA
ACTTTCTCTCTGGTGTATATGGAACCTGTATTGGTGTACAAAATAAATTTGGAGCAGAAGAGAAAAGC
CTTATTGGACTTTCTGGCATTTCATCGCATTGGAGAAATTTTAGGTGGAAGCCTTCTCGGCCTGCTG
AGCAAGAACAATCGTTTTGGTAGAAATCCAGTTGTGCTGTTGGGCATCCTGGTGCATTCATAGCTTTT
TATCTAATATTTCTCAACATGCCTGGAGATGCCCGATTGCTCCTGTTAAAGGAAGTACAGCAGTGTCT
TACATCAAATCCAGCAAAGAAGTTGCCATTCTCTGCAGTTTTCTGTTGGGCTTGGAGACAGCTGCTTT
AATACCCAGCTGCTTAGTATCTTGGGCTTCTGTATTCTGAAGACAGCGCCCAGCATTTGCCATCTTC
AAGTTTGTTCAGTCTATTTGCGCAGCCGTGGCATTTTTCTACAGCAACTACCTTCTCCTTACTGGCAA
CTCCTGGTCATGGTGATATTTGGGTTTTTTGGAACAATTTCTTTCTTCACTGTGGAATGGGAAGCTGCC
GCCTTTGTAGCCCGGGCTCTGACTACCGAAGTATCTGA
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:


ACCN: NM_001242533

Insert Size: 1350 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001242533.1](#)

RefSeq Size: 2988 bp

RefSeq ORF: 1350 bp

Locus ID: 79157

UniProt ID: [O43934](#)
Cytogenetics: 17q25.1-q25.2
Protein Families: Transmembrane
MW: 49.2 kDa