

Product datasheet for **SC329984**

MFSD11 (NM_001242536) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: MFSD11 (NM_001242536) Human Untagged Clone
Tag: Tag Free
Symbol: MFSD11
Synonyms: ET
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC329984 representing NM_001242536.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

ATGTCCCCGGAATCTAAAAAGCTTTTCAACATCATTATTTTAGGAGTTGCCTTTATGTTTATGTTCACT
GCCTTTCAAACCTGTGGAAATGTGGCGCAAACGTGCATCAGGAGCTTAAATAGGACAGATTTTCACGGC
AGTGGATATACCAGCATGGCTATTATCTATGGAGTGTCTCTGCTTCAAATTTGATTACACCGTCAGTG
GTTGCCATTGTAGGACCTCAACTCTCTATGTTTGCCAGTGGTTTATTTTACAGCATGTACATTGCCGTT
TTCATCCAGCCTTCCCCTGGTCTTCTACACAGCCTCTGTTTTCATTGGAATTGCTGCTGCTGAGAGT
GACCGAAGAACAGTGTTTATTGCCCTAACGGTGATTAGCCTTGTTGGGACAGTTCTATTCTTCTCATT
CGGAAACCAGATTCTGAAAATGTCTAGGAGAAGATGAGTCTTCTGATGACCAGGACATGGAAGTCAAC
GAGTCTGCCCAGAACAATCTGACAAAGGAGTAGATGCTTTTAAAAAGTCTTTTAAAGTTATGTGTCACC
AAGGAGATGCTCCTTCTTAGTATTACAACCTGCTTATACAGGTCTGGAATTAACCTTCTTCTCTGGTGTA
TATGGAACCTGTATTGGTGCTACAAATAAATTTGGAGCAGAAGAGAAAAGCCTTATTGGACTTTCTGGC
ATTTTCATCGGCATTGGAGAAAATTTAGGTGGAAGCCTTTCGGCCTGCTGAGCAAGAACAATCGTTTT
GGTAGAAATCCAGTTGTGCTGTTGGCATCCTGGTGCACCTCATAGCTTTTTATCTAATATTTCTCAAC
ATGCCTGGAGATGCCCGATTGCTCCTGTTAAAGGAACGTGACAGCAGTGTACATCAAATCCAGCAAA
GAAGTTGCCATTCTGTCAGTTTTCTGTTGGCCTTGGAGACAGCTGCTTAATACCCAGCTGCTTAGT
ATCTTGGGCTTTCTGTATTCTGAAGACAGCGCCCGCAGCATTGGCCATCTTCAAGTTTGTTAGTCTATT
TGCGCAGCCGTGGCATTCTTCTACAGCAACTACCTTCTCCTTCACTGGCAACTCCTGGTCTAGGTGATA
TTTGGGTTTTTTGGAACAATTTCTTCTTCACTGTGGAATGGGAAGCTGCCGCTTTGTAGCCCGCGC
TCTGACTACCGAAGTATCTGA
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001242536

Insert Size: 1194 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001242536.1](#)

RefSeq Size: 2832 bp

RefSeq ORF: 1194 bp

Locus ID: 79157
UniProt ID: [O43934](#)
Cytogenetics: 17q25.1-q25.2
Protein Families: Transmembrane
MW: 43.3 kDa