

Product datasheet for **SC330000**

AHCYL1 (NM_001242674) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: AHCYL1 (NM_001242674) Human Untagged Clone
Tag: Tag Free
Symbol: AHCYL1
Synonyms: DCAL; IRBIT; PPP1R78; PRO0233; XPVKONA
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC330000 representing NM_001242674.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

ATGCAGGAGTTCACCAAATTCACCAAAAACCTGGCCGAAGATCTTTGTCTCGCTCGATCTCACAGTCC
TCCACTGACAGCTACAGTTCAGCTGCATCTACACAGATAGCTCTGATGATGAGGTTTCTCCCGAGAG
AAGCAGCAAACCAACTCCAAGGGCAGCAGCAATTTCTGTGTGAAGAACATCAAGCAGGCAGAATTTGGA
CGCCGGGAGATTGAGATTGCAGAGCAAGACATGTCTGCTCTGATTTCACTCAGGAAACGTGCTCAGGGG
GAGAAGCCCTTGGCTGGTGTCTAAAATAGTGGGCTGTACACACATCACAGCCAGACAGCGGTGTTGATT
GAGACACTCTGTGCCCTGGGGCTCAGTCCGCTGGTCTGCTTGAACATCTACTCAACTCAGAATGAA
GTAGCTGCAGCACTGGCTGAGGCTGGAGTTGCAGTGTTCGCTTGAAGGGCGAGTCAGAAGTACTTC
TGGTGGTGTATTGACCGCTGTGTGAACATGGATGGGTGGCAGGCCAACATGATCCTGGATGATGGGGGA
GACTTAACCCACTGGGTTTATAAGAAGTATCCAAACGTGTTTAAAGAAGATCCGAGGCATTGTGGAAGAG
AGCGTGACTGGTGTTCACAGGCTGTATCAGCTCTCAAAGCTGGGAAGCTCTGTGTTCCGGCCATGAAC
GTCAATGATTCTGTTACCAAACAGAAGTTTGATAAATTGACTGCTGCCGAGAATCCATTTGGATGGC
CTGAAGAGGACCACAGATGTGATGTTTGGTGGGAAACAAGTGGTGGTGTGTTGGCTATGGTGAGGTAGGC
AAGGGCTGCTGTGCTGCTCTCAAAGCTCTGGAGCAATTGTCTACATTACCGAAATCGACCCCATCTGT
GCTCTGCAGGCCCTGCATGGATGGGTTTCAAGGTTGTAAGCTAAATGAAGTCATCCGGCAAGTCGATGTC
GTAATAAATTGCACAGGAAATAAGAATGTAGTGACACGGGAGCACTTGGATCGCATGAAAAACAGTTGT
ATCGTATGCAATATGGGCCACTCCAACACAGAAATCGATGTGACCAGCCTCCGCACTCCGGAGCTGACG
TGGGAGCGAGTACGTTCTCAGGTGGACCATGTCATCTGGCCAGATGGCAAACGAGTTGCTCCTCTGGCA
GAGGGTCGTCTACTCAATTTGAGCTGCTCCACAGTTCCACCTTTGTTCTGTCCATCACAGCCACAACA
CAGGCTTTGGCACTGATAGAATCTATAATGCACCCGAGGGGCGATACAAGCAGGATGTGTAATTGCTT
CCTAAGAAAATGGATGAATACGTTGCCAGCTTGCATCTGCCATCATTGATGCCACCTTACAGAGCTG
ACAGATGACCAAGCAAAATATCTGGGACTCAACAAAAATGGGCCATTCAAACCTAATTATTACAGATAC
TAA
  
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:


ACCN: NM_001242674

Insert Size: 1452 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001242674.1](#)

RefSeq Size: 3690 bp

RefSeq ORF: 1452 bp

Locus ID: 10768

UniProt ID: [O43865](#)

Cytogenetics:	1p13.3
Protein Families:	Druggable Genome
Protein Pathways:	Cysteine and methionine metabolism, Metabolic pathways, Selenoamino acid metabolism
MW:	53.8 kDa
Gene Summary:	<p>The protein encoded by this gene interacts with inositol 1,4,5-trisphosphate receptor, type 1 and may be involved in the conversion of S-adenosyl-L-homocysteine to L-homocysteine and adenosine. Several transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Jun 2011]</p> <p>Transcript Variant: This variant (3) differs in the 5' UTR and coding sequence compared to variant 1. The resulting isoform (b) is shorter at the N-terminus compared to isoform a. Variants 2, 3, 4, and 5 all encode the same isoform (b). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>