

Product datasheet for SC330008

TMEM139 (NM_001242774) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	TMEM139 (NM_001242774) Human Untagged Clone
Tag:	Tag Free
Symbol:	TMEM139
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC330008 representing NM_001242774. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGTGCCAATGCACTTACTGGGGAGACTGGAGAAGCCGTTCTCCTCCTGTGCTGCGCCTCCTTCTCA
CTGGGGCTGGCTTTGCTGGGCATAAAGACGGACATCACCCCGTTGCTTATTTCTTTCTCACATTGGGT
GGCTTCTTCTTGTTCCTATCTCCTGGTCCGTTTCTGGAATGGGGGCTTCGGTCCCAGCTCCAATCA
ATGCAGACTGAGAGCCCAGGGCCCTCAGGCAATGCACGGGACAATGAAGCCTTTGAAGTGCCAGTCTAT
GAAGAGGCCGTGGTGGGACTAGAATCCCAGTGCCGCCCAAGAGTTGGACCAACCACCCCTACAGC
ACTGTTGTGATACCCCCAGCACCTGAGGAGGAACAACCTAGCCATCCAGAGGGGTCCAGGAGAGCCAAA
CTGGAACAGAGGCGAATGGCCTCAGAGGGTCCATGGCCAGGAAGGAAGCCCTGGAAGAGCTCCAATC
AACCTTCGGCTTCGGGGACCACGGGCTGTGCCACTGCTCCTGATCTGCAGAGCTTGCGGCAGTCCCC
ACATTAGAGCCTCTGACTCCACCCCTGCCTATGATGTCTGCTTTGGTACCCTGATGATGATAGTGT
TTTTATGAGGACAACCTGGGCACCCCTTAA
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



[View online »](#)

ACCN:	NM_001242774
Insert Size:	651 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001242774.2</u>
RefSeq Size:	2189 bp
RefSeq ORF:	651 bp
Locus ID:	135932
UniProt ID:	<u>Q8IV31</u>
Cytogenetics:	7q34
Protein Families:	Transmembrane
MW:	23.7 kDa
Gene Summary:	May be involved in cellular trafficking of proteins such as SLC4A1.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (3) differs in the 5' UTR compared to variant 1. Variants 1, 2, 3, 4, 8, and 9 encode the same protein.