

Product datasheet for **SC330033**

PRKACB (NM_001242860) Human Untagged Clone

Product data:


Product Type: Expression Plasmids
Product Name: PRKACB (NM_001242860) Human Untagged Clone
Tag: Tag Free
Symbol: PRKACB
Synonyms: CAFD2; PKA C-beta; PKACB
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC330033 representing NM_001242860.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGAGTGCACGCAAATCATCAGATGCATCTGCTTCTCCTCTTCAGAAATATCTGATTCCTTTGTGAAA  
GAGTTTCTAGCCAAAGCCAAAGAAGACTTTTTGAAAAAATGGGAGAATCCAACCTCAGAATAATGCCGGA  
CTTGAAGATTTTAAAGGAAAAAAACCCTTGGACAGGTTCAATTTGGAAGAGTCATGTTGGTAAACAC  
AAAGCCACTGAACAGTATTATGCCATGAAGATCTTAGATAAGCAGAAGGTTGTTAACTGAAGCAAATA  
GAGCATACTTTGAATGAGAAAAGAATATTACAGGCAGTGAATTTTCCTTTCCTTGTTCGACTGGAGTAT  
GCTTTTAAAGGATAATTCTAATTTATACATGGTTATGGAATATGTCCTGGGGGTGAAATGTTTTACAT  
CTAAGAAGAATTGGAAGGTTTCAGTGAGCCCATGCACGGTTCTATGCAGCTCAGATAGTGCTAACATTC  
GAGTACCTCCATTCAGTACCTCATCTACAGAGATCTAAAACCTGAAAACTCTTAATTGACCATCAA  
GGCTATATCCAGGTCACAGACTTTGGGTTTGCCAAAAGAGTTAAAGGCAGAACTTGGACATTATGTGGA  
ACTCCAGAGTATTTGGTCCAGAAATAATTCTCAGCAAGGGCTACAATAAGGCAGTGGATTGGTGGGCA  
TTAGGAGTGCTAATCTATGAAATGGCAGCTGGCTATCCCCATTCTTTCAGACCAACCAATTGAGATT  
TATGAAAAGATTGTTTCTGAAAGGTCCGATCCCATCCCACTTCAGTTCAGATCTCAAGGACCTTCTA  
CGAACCTGCTGCAGGTGGATTTGACCAAGAGATTTGAAATCTAAAGAATGGTGTGAGTATATAAAA  
ACTCACAAGTGGTTTCCACGACAGATTGGATTGCTATTTACCAGAGGAAGTTGAAGCTCCATTCCATA  
CCAAAGTTTAGAGGCTCTGGAGATACCAGCAACTTTGATGACTATGAAGAAGAAGATATCCGTGTCTCT  
ATAACAGAAAAATGTGCAAAAAGAATTTGGTGAATTTAA
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:


ACCN: NM_001242860

Insert Size: 1074 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001242860.1](#)

RefSeq Size: 4324 bp

RefSeq ORF:	1074 bp
Locus ID:	5567
UniProt ID:	P22694
Cytogenetics:	1p31.1
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	Apoptosis, Calcium signaling pathway, Chemokine signaling pathway, Dilated cardiomyopathy, Gap junction, GnRH signaling pathway, Hedgehog signaling pathway, Insulin signaling pathway, Long-term potentiation, MAPK signaling pathway, Melanogenesis, Olfactory transduction, Oocyte meiosis, Prion diseases, Progesterone-mediated oocyte maturation, Taste transduction, Vascular smooth muscle contraction, Vibrio cholerae infection, Wnt signaling pathway
MW:	41.3 kDa
Gene Summary:	<p>The protein encoded by this gene is a member of the serine/threonine protein kinase family. The encoded protein is a catalytic subunit of cAMP (cyclic AMP)-dependent protein kinase, which mediates signalling through cAMP. cAMP signaling is important to a number of processes, including cell proliferation and differentiation. Multiple alternatively spliced transcript variants encoding distinct isoforms have been observed. [provided by RefSeq, Jul 2014]</p> <p>Transcript Variant: This variant (7) differs in the 5' exon structure, compared to variant 1. The encoded isoform (7, also known as Cbeta4abc) has a distinct N-terminus and is shorter than isoform 1.</p>