

## Product datasheet for **SC330059**

### **ZA20D3 (ZFAND6) (NM\_001242911) Human Untagged Clone**

#### Product data:


**Product Type:** Expression Plasmids  
**Product Name:** ZA20D3 (ZFAND6) (NM\_001242911) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** ZFAND6  
**Synonyms:** AWP1; ZA20D3; ZFAND5B  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC330059 representing NM\_001242911.  
**Blue**=Insert sequence **Red**=Cloning site **Green**=Tag(s)

```
ATGGCTCAAGAACTAATCACAGCCAAGTGCCTATGCTTTGTTCCACTGGCTGTGGATTTTATGAAAC
CCTCGTACAAATGGCATGTGTTTCAGTATGCTATAAAGAACATCTTCAAAGACAGAATAGTAGTAATGGT
AGAATAAGCCACCTGCAACCTCTGTCTAGTCTGTCTGAATCTTTACCAGTTCAATGCACAGATGGC
AGTGTGCCAGAAGCCAGTCAGCATTAGACTCTACATCTTCATCTATGCAGCCAGCCCTGTATCAAAT
CAGTCACTTTTATCAGAATCTGTAGCATCTTCTCAATTGGACAGTACATCTGTGGACAAAGCAGTACCT
GAAACAGAAGATGTGCAGGCTTCAGTATCAGACACAGCACAGCAGCCATCTGAAGAGCAAAGCAAGTCT
CTTGAAAAACCGAAACAAAAAAGAATCGCTGTTTCATGTGCAGGAAGAAAGTGGGACTTACTGGGTTT
GAATGCCGGTGTGAAATGTTTACTGTGGTGTACACCGTTACTCAGATGTACACAATTGCTTTACAAT
TACAAAGCCGATGCTGCTGAGAAAATCAGAAAAGAAAATCCAGTAGTTGTTGGTAAAAAGATCCAAAAG
ATTTGA
```

**Restriction Sites:** Sgfl-Mlul



[View online »](#)

**Plasmid Map:**


**ACCN:** NM\_001242911

**Insert Size:** 627 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001242911.1](#)

**RefSeq Size:** 1664 bp

**RefSeq ORF:** 627 bp

**Locus ID:** 54469

**UniProt ID:** [Q6FIF0](#)

**Cytogenetics:** 15q25.1

**MW:** 22.6 kDa

**Gene Summary:** Involved in regulation of TNF-alpha induced NF-kappa-B activation and apoptosis. Involved in modulation of 'Lys-48'-linked polyubiquitination status of TRAF2 and decreases association of TRAF2 with RIPK1. Required for PTS1 target sequence-dependent protein import into peroxisomes and PEX5 stability; may cooperate with PEX6. In vitro involved in PEX5 export from the cytosol to peroxisomes (By similarity).[UniProtKB/Swiss-Prot Function]  
Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1-7 all encode isoform a.