

## Product datasheet for **SC330068**

### ZNF410 (NM\_001242924) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** ZNF410 (NM\_001242924) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** ZNF410  
**Synonyms:** APA-1; APA1  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC330068 representing NM\_001242924.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGTTATCAGATGAGTTAGAATCCAAACCAGAGCTCCTGGTACAGTTTGTTCAGAATACGTCCATCCCA
TTGGGACAGGGGCTTGTAGAATCAGAAGCTAAAGATATTACTTGCTTGTCCCTCCTTCCCGTGACTGAA
GCCTCAGAATGCAGTCGGCTAATGTTACCAGAACCTGTTCCATGGAGAGAAGAGGATGGAAGTCAAGC
TGCAGTGACTTGAATGATACTACAAATCATTCTAACTCCTCAAGGAGGTCCTCCTCAGCTGTTTTG
AGAAGCCTTCGGGTGAATGTGGGTCCAGACGGAGAGGAGACGAGAGCTCAGACTGTACAGAAATCCCGG
GAGTTTTTGTCCACTTCAGAGTCTTCTAGCTTGTGCAAGATCTACAGCCAAGTGATAGCACTTCTTTT
ATTCTTCTAACCTAACAAAGAGCAGGTCTGGGCTCTTCAGCTGAGCACTTAGTGTTGTACAGGATGAG
GCAGAAGATTCAGGGAATGATTTCTCCTCCAGTGAGAGCACAGACAGTAGCATTCCATGGTTCTCCGG
GTTCAGGAGTTGGCCATGACAGTTTGATTGCTGCTACTCGTGCACAACCTGGCAAAGAATGCAAAAACC
AGCAGCAATGGAGAAAATGTCCACCTTGGTTCTGGTGATGGGCAGTCAAAAGATTCTGGGCCCTTCTCT
CAAGTGGAAAAGAAGCTCAAGTGTACAGTTGAAGGTTGTGACCGGACATTTGTATGGCCAGCTCACTTT
AAATACCACCTCAAGACTCATCGAAATGACCGCTCCTTCATCTGTCTGCAGAAGGTTGTGGGAAAAGC
TTCTATGTGCTGCAGAGGCTGAAGTGCACATGAGGACCCACAATGGAGAGAAGCCCTTTATGTCCAT
GAGTCTGGCTGTGGTAAGCAGTTTACTACAGCTGGAACCTGAAGAACCACCGGCGCATCCACACAGGA
GAGAAACCTTTCCTTTGTGAAGCCCAAGGATGTGGCCGTTTCCTTTGCTGAGTATTCTAGCCTCCGAAA
CATCTGGTGGTTCACTCAGGAGAGAAGCCTCATCAGTGCCAAGTCTGTGGGAAGACCTTCTCTCAGAGT
GGAAGCAGGAATGTGCATATGAGAAAGCATCACCTGCAGCTGGGAGCAGCTGGGAGTCAAGAGCAGGAG
CAAAGTCTGAGCCACTAATGGGAGTAGTTTGTCTGAAGAGGCTTCAAGTACCCAGTAAAAACCTGGTG
TCTATGAATCCAGCCAGCCTTGGTGGAGAGTCTTGAACCTACCAAATACCAATTCTATCCTGGGA
GTTGATGATGAGGTGCTTGTGAAGGATCCCCACGTTCCCTGTCTTCAAGTGCCTGATGTGACACATCAC
CTGGTGACCATGCAGTCAGGAGGCAATCATATGAAGTTTCTGTCTTAACTGCAGTAAATCCACAAGAG
AGTCTCGCTCCGTTGCCAGGCTGGAGTGCAGTGGTGCCTTCTCGGCTCACTGCAACCTCTGCCTCCCA
GGTTCAAGCGATTCTCCTGCCTCAGCCTCTGA
  
```

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001242924  
**Insert Size:** 1551 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001242924.1</u>
<b>RefSeq Size:</b>	2814 bp
<b>RefSeq ORF:</b>	1551 bp
<b>Locus ID:</b>	57862
<b>UniProt ID:</b>	<u>Q86VK4</u>
<b>Cytogenetics:</b>	14q24.3
<b>Protein Families:</b>	Transcription Factors
<b>MW:</b>	55.8 kDa
<b>Gene Summary:</b>	Transcription factor that activates transcription of matrix-remodeling genes such as MMP1 during fibroblast senescence.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (a).