

Product datasheet for **SC330143**

PDIK1L (NM_001243533) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PDIK1L (NM_001243533) Human Untagged Clone
Tag:	Tag Free
Symbol:	PDIK1L
Synonyms:	CLIK1L; STK35L2
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC330143 representing NM_001243533. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
ATGGTGAGTAGCCAGCCAAAGTACGATCTAATACGGGAGGTAGGCCGAGGTAGTTACGGTGTGTGTAT
GAAGCAGTCATCAGAAAGACCTCTGCACGGGTGGCAGTGAAGAAAATTCGATGTCACGCACCTGAAAAT
GTTGAACTAGCCCTTCGTGAGTTCTGGGCACTAAGCAGTATCAAGAGCCAACATCCAATGTGATTCCAC
TTGGAGGAATGCATCCTACAAAAGGATGGGATGGTGCAAAGATGTCCCACGGCTCTAATTCTTCCCTT
TATTTACAGCTTGTAGAACTTCATTAAGGAGAAAATGCCTTTGATCCCAGAAGCGCCTATTATTTG
TGGTTTGTGATGGATTTTGTGACGGAGGAGATATGAATGAGTATCTGTTGTCCAGGAAACCAATCGT
AAAACAAACACCAGCTTCATGCTTCAGCTGAGCAGTGCCCTGGCTTTCTTGCATAAAAACAGATCATC
CACCGAGATCTTAAGCCTGATAACATCCTGATTTCTCAAACAGGTTGGATACCAGTGACTTGGAACCT
ACCTCAAAGTGGCTGATTTTGGTCTAAGTAAAGTTTGTTCAGCCTCTGGGCAGAACCCAGAAGAACCT
GTCAGTGTAAACAAGTGTTCCTTTCCACAGCATGTGGAACAGATTTTACATGGCTCCTGAAGTTTGG
GAAGGACATTACACAGCAAAAGCTGACATCTTTGCTCTGGGGATTATCATCTGGGCAATGCTGGAAGG
ATCACATTCATAGACACAGACAAAGAAGGAACCTTTGGGGAGTTATGTAAAACAAGGAACCTGAGATT
GTGCCTGTTGGGGAGGCACTTCTGAAAAATCCAAAATGGAACCTCTCATTCTGTGAAGAAAAATCT
ATGAATGGGCGAATGAAACAACGATTAAGGAAATGCTGGCTGCAAACCCCTCAGGATCGTCCAGATGCT
TTTGAAC TAGAACTCAGATTAGTACAAAATTGCATTTAAAGATAGCAGCTGGGAAACGTGA
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:



ACCN: NM_001243533

Insert Size: 1026 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001243533.1](https://www.ncbi.nlm.nih.gov/nuccore/NM_001243533.1)

RefSeq Size: 4344 bp
RefSeq ORF: 1026 bp
Locus ID: 149420
UniProt ID: [Q8N165](#)
Cytogenetics: 1p36.11
Protein Families: Druggable Genome, Protein Kinase
MW: 38.5 kDa