

Product datasheet for **SC330184**

PAMM (FAM213A) (NM_001243780) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: PAMM (FAM213A) (NM_001243780) Human Untagged Clone
Tag: Tag Free
Symbol: PRXL2A
Synonyms: Adrx; C10orf58; FAM213A; PAMM
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC330184 representing NM_001243780.
Blue=Insert sequence **Red**=Cloning site **Green**=Tag(s)


```
ATGTCTTTCTCCAGGACCCAAGTTTCTTCACCATGGGGATGTGGTCCATTGGTGCAGGAGCCCTGGGG
GCTGCTGCCTTGGCATTGCTGCTTGCCAACACAGACGTGTTTCTGTCCAAGCCCCAGAAAGCGGCCCTG
GAGTACCTGGAGGATATAGACCTGAAAACACTGGAGAAGGAACCAAGGACTTTCAAAGCAAAGGAGCTA
TGGGAAAAAATGGAGCTGTGATTATGGCCGTGCGGAGGCCAGGCTGTTTCTCTGTGCGAGAGGAAGCT
GCGGATCTGTCTCCCTGAAAAGCATGTTGGACCAGCTGGGCGTCCCCCTCTATGCAGTGGTAAAGGAG
CACATCAGGACTGAAGTGAAGGATTTCCAGCCTTATTTCAAAGGAGAAATCTTCTGGATGAAAAGAAA
AAGTTCTATGGTCCACAAGGCGGAAGATGATGTTTATGGGATTTATCCGTCTGGGAGTGTGGTACAAAC
TTCTTCCGAGCCTGGAACGGAGGCTTCTCTGGAACCTGGAAGGAGAAGGCTTCATCCTTGGGGGAGTT
TTCGTGGTGGGATCAGGAAAGCAGGGCATTCTTCTTGGACCCGAGAAAAAGAATTTGGAGACAAAGTA
AACCTACTTTCTGTTCTGGAAGCTGCTAAGATGATCAAACCACAGACTTTGGCCTCAGAGAAAAATGA
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:



ACCN: NM_001243780

Insert Size: 690 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001243780.1
RefSeq Size:	2239 bp
RefSeq ORF:	690 bp
Locus ID:	84293
UniProt ID:	Q9BRX8
Cytogenetics:	10q23.1
Protein Families:	Secreted Protein, Transmembrane
MW:	25.8 kDa
Gene Summary:	<p>Involved in redox regulation of the cell (PubMed:26438880, PubMed:19951071). Acts as an antioxidant (PubMed:19951071, PubMed:26438880). Inhibits TNFSF11-induced NFKB1 and JUN activation and osteoclast differentiation (PubMed:19951071). May affect bone resorption and help to maintain bone mass (PubMed:19951071). Acts as a negative regulator of macrophage-mediated inflammation by inhibiting macrophage production of inflammatory cytokines, probably through suppression of the MAPK signaling pathway (PubMed:26438880). [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (4) uses an alternate splice site in the 5' UTR compared to variant 1. Variants 1, 2, 3 and 4 encode the same isoform (1).</p>