

## Product datasheet for **SC330318**

### AP4S1 (NM\_001254726) Human Untagged Clone


#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** AP4S1 (NM\_001254726) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** AP4S1  
**Synonyms:** AP47B; CLA20; CLAPS4; CPSQ6; SPG52  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC330318 representing NM\_001254726.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

ATGATAAAATTTTCTCATGGTGAATAAACAAGGGCAGACTCGACTTTCTAAGTACTATGAACATGTG  
 GATATTAATAAGCGTACACTTCTGGAAACAGAAGTCATAAAGAGCTGTCTCTCGATCCAATGAACAA  
 TGCTCTTTCATTGAATATAAGGATTTTAAGCTGATATATCGGCAGTATGCAGCTCTTTCATTGTGGTT  
 GGAGTTAATGACACTGAGAACGAGATGGCTATTTATGAATTCATTCATACTTTGTGGAAGTTTTAGAT  
 GAGTATTTAGCCGAGTGAACCAATTGATGAACCTCCAAAATATGCTCAGCCCTGGAGCCCCAGCAG  
 ACTTGCTTTTCTCCAGACAGTTCATCATTCAAAGGAGCTGCCTCCACCACCCCATCTACTGA

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**



**ACCN:** NM\_001254726



[View online >](#)

<b>Insert Size:</b>	408 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001254726.1</a>
<b>RefSeq Size:</b>	1732 bp
<b>RefSeq ORF:</b>	408 bp
<b>Locus ID:</b>	11154
<b>UniProt ID:</b>	<a href="#">Q9Y587</a>
<b>Cytogenetics:</b>	14q12
<b>Protein Pathways:</b>	Lysosome
<b>MW:</b>	15.7 kDa
<b>Gene Summary:</b>	<p>This gene encodes a member of the adaptor complexes small subunit protein family. These proteins are components of the heterotetrameric adaptor protein complexes, which play important roles in the secretory and endocytic pathways by mediating vesicle formation and sorting of integral membrane proteins. The encoded protein is the small subunit of adaptor protein complex-4, which is associated with both clathrin- and nonclathrin-coated vesicles. Mutations in this gene are associated with spastic quadriplegic cerebral palsy-6. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 6. [provided by RefSeq, Dec 2011]</p> <p>Transcript Variant: This variant (3) lacks an exon in the 3' coding region, but maintains the reading frame, compared to variant 1. The encoded isoform (3) is shorter than isoform 1.</p> <p>Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>