

## Product datasheet for **SC330348**

### **PNPLA8 (NM\_001256010) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	PNPLA8 (NM_001256010) Human Untagged Clone
Tag:	Tag Free
Symbol:	PNPLA8
Synonyms:	IPLA2-2; IPLA2G; iPLA2gamma; MMLA; PNPLA-gamma
Vector:	pCMV6-Entry (PS100001)




[View online »](#)

Fully Sequenced ORF: >SC330348 representing NM\_001256010.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGTCCCGTATTAAGTACTTTGAACTCTGTTTCAAAGGCTGTTTTGGCAATCAAAATGAAATGATT
TCACGTTTAGCTCAATTTAAGCCAAGTTCCCAAATTTAAGAAAAGTATCGGATAGTGGCTGGTAAAA
CAGAAAAACATCAAAACAGCCATCAAATCTCTGAAAAATATAGTGACAAATCAGCAGAAAAGATCCT
TTCCAGAAAGAGAAAAGTACATTATAGACAAAAGAAGATATAGGTAACGCAGTCTTTTTCATTAC
ACAAGTTCTATAACCACAAAATTTGGAGACTCATTCTACTTTTTATCAAATCATATTAATTCATATTC
AAACGTAAAGGAAAAATGTCTCAACAAAAGGAAAAATGAACATTTCCGGGACAAATCAGAACTTGAAGAT
AAAAAGGTAGAAGAGGGGAAATTAAGATCTCCAGATCCTGGCATCCTGGCTTATAAGCCAGGCTCAGAA
TCTGTACATACGGTGGACAAGCCTACAAGTCTTCTGCGATACCTGATGTTCTTCAAGTTTCAACTAAA
CAAAGTATTGCTAACTTTCTTCTCGTCCACGGAAGGTGTACAAGCTTTAGTAGGTGGTTATATTGGT
GGACTTGTCCCAAATTAAGTATGATTCAAAGAGTCAGTCAGAAGAACAGGAAGAGCCTGCTAAAAC
GATCAGGCTGTCAGCAAAGACAGAAATGCAGAGGAGAAAAAGCGTTTATCTCTTACGAGAAAAGATT
ATCGCAAGGGTGAGTATTGATAACAGGACCCGGGCATTAGTTCAGGCATTAAGAAGAACAACCTGACCCA
AAGCTCTGCATTACTAGGGTTGAAGAACTGACTTTTCATCTTCTAGAATTCCTGAAGGAAAAGGAGTG
GCTGTCAAGGAAAAGAAATTATCCATATTTATTACGACTGAGACAAATTAAGGATGAAACTCTTCAGGCT
GCAGTTAGAGAAATTTGGCCCTAATTGGCTATGTGGATCCAGTAAAAGGGAGAGGAATCCGAATTCCTC
TCAATTGATGGTGGAGGAACAAGGGCGTGGTTGCTCTCCAGACCCTACGAAAATTAGTTGAACTTACT
CAGAAGCCAGTTCATCAGCTCTTTGATTACATTTGGTGTAAAGCACAGGTGCCATATTAGCTTTCATG
TTGGGGTTGTTTCATATGCCCTTGGATGAATGTGAGGAACTTTATCGAAAATTAGGATCAGATGTATTT
TCACAAAATGTCATTGTTGGAACAGTAAAAATGAGTTGGAGCCATGCATTTTATGACAGTCAAACATGG
GAAAACATTCCTTAAGGATAGGATGGGATCTGCACTGATGATTGAAACAGCAAGAAAACCCACATGTCCT
AAGGTAGCTGCTGTAAGTACCATAGTAAATAGAGGGATAACACCCAAAGCTTTTGTGTTCAAGAACTAT
GGTCATTTTCTGGAATCAACTCTCATTATTTGGGAGGCTGTGAGTATAAAATGTGGCAGGCCATTAGA
GCCTCATCTGCTGCTCCAGGCTACTTTGCAGAATATGCATTGGGAAATGATCTTCATCAAGATGGAGGT
TTGCTTGAATAACCTTCGGCATTAGCTATGCATGAGTGAAATGTCTTTGGCCAGATGTGCCGTTA
GAGTGATAGTATCCCTGGGCACTGGACGTTATGAGAGTGATGTGAGAAACACGGTAACATACACAAGC
TTGAAAACATAAATTTCTAATGTTATCAACAGTGCTACAGATACAGAAGAAGTCCATATAATGCTTGAT
GGCCTGTACCTCCTGACACCTATTTAGATTCATCCTGTAATGTGTGAAAACATACCTCTAGATGAA
AGTCGAAATGAAAAGCTGGATCAGCTGCAGTTGGAAGGGTTGAAATACATAGAAAAGAAATGAACAAAA
ATGAAAAAAGTTGAAAAATATTAAGTCAAGAAAAACAACCTCTGCAGAAAATTAATGATTGGATAAAA
TTAAAACTGATATGATGAAGGACTTCCATTCTTTTCAAATTTGA
```

Restriction Sites: SgfI-MluI

## Plasmid Map:



ACCN: NM\_001256010

Insert Size: 2049 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001256010.1](#)

RefSeq Size: 4574 bp

RefSeq ORF: 2049 bp

Locus ID: 50640

UniProt ID: [Q9NP80](#)

Cytogenetics: 7q31.1

MW: 77 kDa

**Gene Summary:**

This gene encodes a member of the patatin-like phospholipase domain containing protein family. Members of this family are phospholipases which catalyze the cleavage of fatty acids from membrane phospholipids. The product of this gene is a calcium-independent phospholipase. Mutations in this gene have been associated with mitochondrial myopathy with lactic acidosis. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2015]

Transcript Variant: This variant (5) differs in the 5' UTR and coding region and uses a downstream start codon compared to variant 1. The resulting protein (isoform 3) has a shorter N-terminus compared to isoform 1. Variants 5 and 6 encode the same protein.

Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.