

Product datasheet for **SC330355**

Contactin 1 (CNTN1) (NM_001256063) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Contactin 1 (CNTN1) (NM_001256063) Human Untagged Clone
Tag:	Tag Free
Symbol:	Contactin 1
Synonyms:	F3; GP135; MYPCN
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC330355 representing NM_001256063. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGAAATGTGGTTGCTGGTCAGTCATCTTGTGATAATATCTATTACTACCTGTTTAGCAGAGTTTACA
TGGTATAGAAGATATGGTCATGGAGTTTCTGAGGAAGACAAAGGATTTGGACCAATTTTTGAAGAGCAG
CCAATCAATACCATTTATCCAGAGGAATCACTGGAAGGAAAAGTCTCACTCAACTGTAGGGCAGGAGCC
AGCCCTTCCCGGTTTACAAATGGAGAATGAATAATGGGGACGTTGATCTCACAAGTGATCGATACAGT
ATGGTAGGAGGAAACCTTGTATCAACAACCTGACAAACAGAAAGATGCTGGAATATACTACTGTTTA
GCATCTAATAACTACGGGATGGTCAGAAGCACTGAAGCAACCCTGAGCTTTGGATATCTTGATCCTTTC
CCACCTGAGGAACGTCCTGAGGTCAGAGTAAAAGAAGGAAAGGAATGGTGCTTCTCTGTGACCCCCCA
TACCATTTCCAGATGATCTTAGCTATCGCTGGCTTCTAAATGAATTTCCGTGATTTATCACAATGGAT
AAACGGCGATTTGTGCTCAGACAAATGGCAATCTCTACATTGCAAAATGTTGAGGCTTCCGACAAAGGC
AATTATTCCTGCTTTGTTCCAGTCCTTCTATTACAAAGAGCGTGTTTCAGCAAATTCATCCCCTCATT
CCAATACCTGAACGAACAACAAACCATATCCTGCTGATATTGTAGTTCAGTTCAAGGATGTATATGCA
TTGATGGGCCAAAATGTGACCTTAGAATGTTTTGCACTTGAAATCCTGTTCCGGATATCCGATGGCGG
AAGTTCTAGAACCAATGCCAAGCACTGCTGAGATTAGCACCTCTGGGGCTGTTCTTAAGATCTTCAAT
ATTCAGCTAGAAGATGAAGGCATCTATGAATGTGAGGCTGAGAACATTAGAGGAAAGGATAAACATCAA
GCAAGAATTTATGTTCAAGCATTCCCTGAGTGGGTAGAACACATCAATGACACAGAGGTGGACATAGGC
AGTGATCTCTACTGGCCTTGTGTGGCCACAGGAAAGCCATCCCTACAATCCGATGGTTGAAAAATGGA
TATGCGTATCATAAAGGGGAATTAAGACTGTATGATGTGACTTTTGAAAAAGCCGGAATGTATCAGTGC
ATAGCTGAAAACACATATGGAGCCATTTATGCAAAATGCTGAGTTGAAGATCTTGGCGTTGGCTCCAAT
TTTGAAATGAATCCTATGAAGAAAAAGATCCTGGCTGCTAAAGGTGGAAGGGTGATAATTGAATGCAAA
CCTAAAGCTGCACCGAAACCAAGTTTTTATGGAGTAAGGGACAGAGTGGCTTGTCAATAGCAGCAGA
ATACTCATTTGGGAAGATGGTAGCTTGGAAATCAACAACATTACAAGGAATGATGGAGGTATCTATACA
TGCTTTGCAGAAAATAACAGAGGGAAGCTAATAGCACTGGAACCCTTGTATCACAGATCTACGCGA
ATTATATTGGCCCAATTAATGCCGATATCACAGTTGGAGAAAACGCCACCATGCAGTGTGCTGCGTCC
TTTGATCCTGCTTGGATCTCACATTTGTTGGTCTTCAATGGCTATGTGATCGATTTTAAACAAAGAG
AATATTCACTACCAGAGGAATTTATGCTGGATTCCAATGGGGAATTAATAATCCGAAATGCGCAGCTG
AAACATGCTGGAAGATACACATGCACTGCCAGACAATTGTGGACAATTTTCAGCTTTCAGCTGACCTT
GTAGTGAGAGGTAACAGAAAAGGTGGAGAGAAAACATGGTGGATTCTTTCTCCCGTCTGTGCA
AGCCTGCCTCAACCTGGTGA

```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001256063
Insert Size:	1884 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001256063.1
RefSeq Size:	3315 bp
RefSeq ORF:	1884 bp
Locus ID:	1272
UniProt ID:	Q12860
Cytogenetics:	12q12
Protein Pathways:	Cell adhesion molecules (CAMs)
MW:	70.6 kDa
Gene Summary:	<p>The protein encoded by this gene is a member of the immunoglobulin superfamily. It is a glycosylphosphatidylinositol (GPI)-anchored neuronal membrane protein that functions as a cell adhesion molecule. It may play a role in the formation of axon connections in the developing nervous system. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2011]</p> <p>Transcript Variant: This variant (3) lacks several 3' exons and includes an alternate 3' terminal exon, compared to variant 1. The resulting isoform (3) is shorter and has a distinct C-terminus, compared to isoform 1. Variants 3 and 4 encode the same isoform 3. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>