

Product datasheet for SC330375

FSBP (NM_001256141) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FSBP (NM_001256141) Human Untagged Clone
Tag:	Tag Free
Symbol:	FSBP
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC330375 representing NM_001256141. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGTAGGAAAGGCTAGATCTTCCAATTTACCTTATCCGAAAAGCTTGATTTGCTAAAGCTTGTGAAG
CCATATGTGAAAATTCTCGAAGAACACACTAATAAACATTCAGTAATAGTAGAAAAGAATAGATGTTGG
GATATCATAGCAGTTAACTATAATGCAATTGGAGTAGACCGCCCTCCTCGAACAGCACAGGGCCTACGC
ACCTTTTATAAAAGGCTCAAAGAATATGCCAAACAGGAGCTATTGCAGCAAAAAGAGACCCAATCAGAT
TTTAAAAGCAATATTTCTGAGCCAACCAAGAAAGTTATGGAGATGATTCCTCCAGATTTCCAGTTTTTGC
CTGGTAAGAGACAGGAACACATACAAAGTGCAAACCTTGGATGAGGAGGCACAGGCTGGTACCAGTTCA
CTACAGGTAATGTTGGATCACCATCCTGTTGCTATTACAGTGGAGGTGAAGCAAGAAGAAGACATTAAA
CCACCTCTCCACTGGTTTTAAATTCTCAACAGAGTGATACTTTAGAGCAAAGAGAAGAATGAATTA
GTACATGTTATGGAAAGATCCTTGTCGCCTTCACTTTCTCTGTTGATATGAGAATGACATCGTCTCCA
TCTTCTATTCCAAGGAGAGATGATTTTTTTCGGCATGAGAGTGGTGAACACTTTAGGTCACTATTAGGG
TATGATCCTCAGATCCTGCAATGTTGAAAGAGGAGCATCAGATAATTTAGAAAATCAAAAAATTTT
GGATTGTATGTTCAAGGAGAAGAGGGATGGATTGAAAAGAAGGCAGCAGCTAGAGGAAGAGCTACTAAGA
GCAAAAATTGAAGTGGAGAAGCTGAAAGCAATTCGCTTACGGCATGATCTACCTGAATATAACAGTCTC
TAA
  
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001256141
Insert Size:	900 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001256141.1](#)

RefSeq Size: 5418 bp

RefSeq ORF: 900 bp

Locus ID: 100861412

UniProt ID: [O95073](#)

Cytogenetics: 8q22.1

MW: 34.8 kDa