

## Product datasheet for **SC330384**

### ADGRE1 (NM\_001256252) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ADGRE1 (NM_001256252) Human Untagged Clone
Tag:	Tag Free
Symbol:	ADGRE1
Synonyms:	EMR1; TM7LN3
Vector:	pCMV6-Entry (PS100001)




[View online »](#)

Fully Sequenced ORF: >SC330384 representing NM\_001256252.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

ATGCGTGGCTTCAACCTGCTCCTCTTCTGGGGATGTTGTGTTATGCACAGCTGGGAAGGGCACATAAGA  
CCCACACGGAAACCAACACAAAGGGTAATAACTGTAGAGACAGTACCTTGTGCCAGCTTATGCCACC  
TGCACCAATACAGTGGACAGTTACTATTGCGCTTGCAAACAAGGCTTCTGTCCAGCAATGGGCAAAAT  
CACTTCAAGGATCCAGGAGTGCAGTGCAAAGATATCAATGAGTGCCTCACCAGCAGCGTCTGCCCTGAG  
CATTCTGACTGTGTCAACTCCATGGGAAGCTACAGTTGCAGCTGTCAAGTTGGATTTCATCTAGAAAC  
TCCACCTGTGAAGACGTGGATGAATGTGCAGATCCAAGAGCTTGCCAGAGCATGCAACTTGAATAAC  
ACTGTTGGAACTACTCTGTTTCTGCAACCCAGGATTTGAATCCAGCAGTGGCCACTTGAGTTTCCAG  
GGTCTCAAAGCATCGTGTGAAGATATTGATGAATGCAGTAAAATGTGCCCATCAATTCAACATGCACC  
AACACTCTGGGAGTACTTTTGCACCTGCCACCCTGGCTTTGCACCAAGCAATGGACAGTTGAATTTT  
ACAGACCAAGGAGTGAATGTAGAGATATTGATGAGTCCCAGCAAGATCCATCAACTGTGGTCTAAT  
TCTATCTGCACCAATGCCCTGGGCTCCTACAGCTGTGGTGCATTGCAGGCTTTCATCCCAATCCAGAA  
GGCTCCCAGAAAGATGGCAACTCAGCTGCCAAAGGGTTCTCTCAAATGTAAGGAAGATGTGATACCC  
GATAATAAGCAGATCCAGCAATGCCAAGAGGGAACCGCAGTGAAACCTGCATATGTCTCCTTTTGTGCA  
CAAATAAATAACATCTTCAGCGTTCTGGACAAAGTGTGTGAAAATAAAACGACCGTAGTTTCTCTGAAG  
AATACAACGTAGAGCTTTGTCCCTGTGCTTAAACAAAATATCCACGTGGACTAAATTCACCAAGGAAGAG  
ACGTCTCCCTGGCCACAGTCTTCTGGAGAGTGTGGAAGCATGACACTGGCATCTTTTGGAAACCC  
TCAGCAAATATCACTCCGGCTGTTCCGGACGGAATACTTAGACATTGAGAGCAAAGTTATCAACAAAGAA  
TGCAGTGAAGAGAATGTGACGTTGGACTTGGTAGCCAAGGGGATAAGATGAAGATCGGGTGTCCACA  
ATTGAGGAATCTGAATCCACAGAGACCACTGGTGTGGCTTTTGTCTCCTTTGTGGCATGGAATCGGTT  
TTAAATGAGCGCTTCTTCAAAGACCACCAAGGCTCCCTTGACCACCTCTGAGATCAAGTGAAGATGAAT  
TCTCGAGTCGTTGGGGGCATAATGACTGGAGAGAAGAAAGACGGCTTCTCAGATCCAATCATCTACT  
CTGGAGAACATTAGCCAAGCAGAAAGTTTGGAGGCCATCTGTGTTTCTGGAGCACTGATGTGAAG  
GGTGAAGATGGACATCCTTTGGCTGTGTGATCCTGGAAGCTTCTGAGACATATACCATCTGCAGCTGT  
AATCAGATGGCAAATCTTGCCGTTATCATGGCGTCTGGGAGCTCACGATGGACTTTTCTGTACATC  
ATTAGCCATGTAGGCATTATCATCTCCTTGGTGTGCTCGTCTTGGCCATCGCCACCTTTCTGCTGTGT  
CGCTCCATCCGAAATCACAAACCTACCTCCACCTGCACCTCTGCGTGTGTCTCCTTTGGCGAAGACT  
CTCTTCTCGCCGTATACACAAGACTGACAACAAGATGGGCTGCGCCATCATCGCGGGCTTCTGCAC  
TACCTTTTCTTGCCTGCTTCTTCTGGATGCTGGTGGAGGCTGTGATACTGTTCTTGATGGTCAGAAAC  
CTGAAGGTGGTGAATTACTTTCAGCTCTCGCAACATCAAGATGCTGCACATCTGTGCCCTTTGGTTATGGG  
CTGCCGATGCTGGTGGTGGTGTCTCTGCCAGTGTGCAGCCACAGGGCTATGGAATGCATAATCGCTGC  
TGGCTGAATACAGAGACAGGGTTCATCTGGAGTTTCTTGGGGCCAGTTTGCACAGTTATAGTGGTAAAGC  
AAATACTACAACAGCCTGGCGAAGTGTGTTCTGAAGGAGGAGCAAGGAGACCTGCGAGATCTGGAATTT  
CCAGGGACGTGTGCAGCTGAGAGGATCAACTCCCTTCTCTGACCTGGACCTTGTGGATCCTGAGGCAG  
AGGCTTTCCAGTGTTAATGCCGAAGTCTCAACGCTAAAAGACACCAGGTTACTGACCTTCAAGGCCTTT  
GCCAGCTTTCATCCTGGGCTGCTCCTGGTGTGGCATTTCAGATTGGACCTGTGGCAGGTGTC  
ATGGCTTACCTGTTCAACATCATCAACAGCCTGCAGGGGCTTTCATCTTCTCATCCACTGTCTGCTC  
AACGGCCAGGTACGAGAAGAAACAAGAGGTGGATCACTGGGAAGACGAAGCCAGCTCCCAGTCCCAG  
ACCTCAAGGATCTTGTCTCCTCCATGCCATCCGCTTCCAAGACGGGTTAA

Restriction Sites: Sgfl-Mlul

## Plasmid Map:



ACCN: NM\_001256252

Insert Size: 2604 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001256252.1](#)

RefSeq Size: 3104 bp

RefSeq ORF: 2604 bp

Locus ID: 2015

UniProt ID: [Q14246](#)

Cytogenetics: 19p13.3-p13.2

Protein Families: Druggable Genome, Transmembrane

MW: 95.9 kDa

**Gene Summary:**

This gene encodes a protein that has a domain resembling seven transmembrane G protein-coupled hormone receptors (7TM receptors) at its C-terminus. The N-terminus of the encoded protein has six EGF-like modules, separated from the transmembrane segments by a serine/threonine-rich domain, a feature reminiscent of mucin-like, single-span, integral membrane glycoproteins with adhesive properties. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2012]

Transcript Variant: This variant (2) lacks an exon in the 5' coding region and has an additional segment in the 3' coding region, but maintains the reading frame, compared to variant 1. The resulting isoform (2) is shorter than isoform 1.