

Product datasheet for **SC330412**

ELOVL1 (NM_001256399) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: ELOVL1 (NM_001256399) Human Untagged Clone
Tag: Tag Free
Symbol: ELOVL1
Synonyms: CGI-88; IKSHD; Ssc1
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC330412 representing NM_001256399.
Blue=Insert sequence **Red**=Cloning site **Green**=Tag(s)


```
ATGGAGGCTGTTGTGAACCTGTACCAAGAGGTGATGAAGCACGCAGATCCCCGGATCCAGGGCTACCCCT
CTGATGGGGTCCCCCTTGCTAATGACCTCCATTCTCCTGACCTACGTGTACTTCGTTCTCTCACTGGG
CCTCGCATCATGGCTAATCGGAAGCCCTCCAGCTCCGTGGCTTCATGATTGTCTACAACCTCTCACTG
GTGGCACTCTCCCTCTACATTGTCTATGAGTTCCTGATGTGGGGCTGGCTGAGCACCTATACCTGGCGC
TGTGACCCTGTGGACTATTCCAACAGCCCTGAGGCACTTAGGATGGTTCGGGTGGCCTGGCTCTTCCTC
TTCTCCAAGTTCATTGAGCTGATGGACACAGTGATCTTTATTCTCCGAAAGAAAGACGGGCAGGTGACC
TTCTACATGTCTCCATCACTCTGTGCTCCCTGGAGCTGGTGGTGGGGGTAAGATTGCCCCGGGA
GGAATGGGCTCTTTCATGCCATGATAAACTCTCCGTGCATGTCATAATGTACCTGTACTACGGATTA
TCTGCCTTTGGCCCTGTGGCACAACCTACCTTTGGTGGAAAAAGCACATGACAGCCATTACAGCTGATC
CAGTTTGTCTGGTCTCACTGCACATCTCCAGTACTACTTTATGTCCAGCTGTAACACAGTACCCA
GTCATTATTCACCTCATCTGGATGTATGGCACCATCTTCTTCATGCTGTTCTCCAACCTCTGGTATCAC
TCTTATACCAAGGGCAAGCGGCTGCCCGTGCACTTCAGCAAAATGGAGCTCCAGGTATTGCCAAGGTC
AAGGCCAACTGA
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:



ACCN: NM_001256399

Insert Size: 840 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq:	NM_001256399.1
RefSeq Size:	1631 bp
RefSeq ORF:	840 bp
Locus ID:	64834
UniProt ID:	Q9BW60
Cytogenetics:	1p34.2
Protein Families:	Transmembrane
MW:	32.7 kDa
Gene Summary:	<p>Catalyzes the first and rate-limiting reaction of the four reactions that constitute the long-chain fatty acids elongation cycle. This endoplasmic reticulum-bound enzymatic process allows the addition of 2 carbons to the chain of long- and very long-chain fatty acids (VLCFAs) per cycle. Condensing enzyme that exhibits activity toward saturated and monounsaturated acyl-CoA substrates, with the highest activity towards C22:0 acyl-CoA. May participate in the production of both saturated and monounsaturated VLCFAs of different chain lengths that are involved in multiple biological processes as precursors of membrane lipids and lipid mediators. Important for saturated C24:0 and monounsaturated C24:1 sphingolipid synthesis. Indirectly inhibits RPE65 via production of VLCFAs.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) uses a different 5' UTR, compared to variant 1. Variants 1 and 2 encode the same protein (isoform 1).</p>