

Product datasheet for **SC330470**

GIDRP88 (R3HCC1L) (NM_001256619) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GIDRP88 (R3HCC1L) (NM_001256619) Human Untagged Clone
Tag:	Tag Free
Symbol:	GIDRP88
Synonyms:	C10orf28; GIDRP86; GIDRP88; PSORT
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC330470 representing NM_001256619.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

ATGCAGCAAGAATCAGAGAGATGCAGAGTTAGAGCCAGAAGGCCTGACATGGCACTTTATGTACCTAAA
 GCTCGTAGGGGTGCAGTACTCCTTAAGACAGGTGATGAAGAAGAAAGCTGTGGTTACCTAACTCTGTG
 GTGAAAGAAAAGCAAAAAGAAAGTTCTCTCTCCCAAAAAGAAAGCTTTAAAGACAAACCGGAGGCTCGA
 AGACTAAATATCAATCCTGATAGAAAGGAGCATAATTGTAGAGAAGAAAAGAAATCTTCAACAAAATTA
 AGAATGGACACATGCCTTCAAAAACAAATAGAGTTTGTCTAAGAGAGGAACCACTGAATCCAAAGAA
 GTATTATCCCAAGGACAACAGCAGGGAGCCCCAAATGCTGGGGTTATAACTAATGCACCTTTGCAGAGA
 CATTTTAAACCAAGAAGGTGGAGTGTTTGAAGTTGAACTACGGATGTGACAGGACATGAGAGGATA
 CTTCTTTACAGGCCTGTTTAAATCAGCAGGCTCAAGTTCCAAGCAAACCATCCAAAATGTGGAA
 TTCTGTGACTTCAGTAGGCATGAACCTGATGGGAAGCATTGAGACAAAGATTGGAAGGCAGAAAT
 GAACTGATACCAAGTTTTGGAGATACTATAGAGTTTCTAGAGTTTTAGTTCTGTCATGAAACCT
 GAGAATATGATTGTACCAATAAACTAAGCTCTGATTCTGAAATTGTACAACAAAGCATGCAACATCA
 GATGGAATATTGAATCCAGCAGCGGAGGCATCACCCTACTTCTGTTCTGGAAGTCCAGATGGTGTC
 TTTGATCAAACCTTGCCTAGATTTTGAAGTTGAGAGTGTAGGTGGTATAGCCAATAGTACAGGTTTCATC
 TTAGATCAAAAAGATACAGATTCATTCTGCAACTATGGGTACATCTCTCTGTCAGAGAGCACAAAT
 GACACTGTTAGTCCAGTAATGATTAGAGAATGTGAGAAGAATGACAGCACTGCTGATGAGTTACATGTA
 AAGCACGAACCTCCTGATACAGCTGTCCTTGCTCATGAAACACATAGAGATAGTGGATTTAAGAATGTA
 GGTGACATTACCAATAAAGCATGTATGATGGACACTACAGGTATGTCCTGTAGTGATCATGTAACGTG
 GATAGCCCTTATGTAGTTGCAGTTAGAATAGCTGATGAGACCTCTATTAATACACGAAGTTTCTCAAAG
 TTTGTAGGAATGAGTGCAGATGCAACCCCTCTTCATGTAGCTAGAAGTGGGAATGACACTGAAGATTTT
 AGCAACCCCTTCTGCTTGCTCAGATATTTATGGTGAGAGTATTTTCATCTCATTTTACAGAGTCAACAGGA
 AAGTTGATAGAGAGCTTGTGAGATTGTGCTTCTCCTTACCTATAAAAAAGATTGCTGGTAGTAATTAT
 AACACTTTTTTGGACTCTGAACTCAGTATGTTAAATGGGACAAAAGTTCTTTTACAGAGTGCCTGGGC
 ATTGACCTGGGTAGTACTGGTGATACAACAGAAGCATTGCACGAAGTAAAGAACTGCCGAAGAGTTCAAA
 ACAGAAGAGCAAGATGACTCAGGGAGTATAGAATTTGGTGATCTTTTCTGATAGGGAATCATCATCT
 ATGGAAACATCCATCGAACCAAAAGCAACTGAACTTCTCACACAGAGGGAATTACTGCCATTGAGGAG
 AGCTGGGAGTCTATGTTTAAAGATGATGGTACTGCCTGGATCCACGCTTCTACAAGAGGGGATTCTG
 ATGCACATAAAGCCCAGAACTGCTCAAAGTTATCAGGGAATACCAAGAGCAGAGAGAGCATCCAG
 GAACCTAGATCTGATTACTACAATCATGAAGTTCTGATATTGACCTCAGTGATTGTGAATCCCACAT
 GTCATTGAAATTTATGACTTTCCCAAGAATTTCTGACTGAAGACCTTCTACGGGTTTTCTGCAGTTAT
 CAAAAGAAAGGATTTGATATTAATGGGTGGATGATACACATGCCTAGGAGTATTCTCCAGTCCAATT
 ACAGCTCGTGATGCGTTGGGTATTAACACACCATGGTGAAGATTCGTCCTTGTACAGGCCACAAGA
 GCAGCCAGGCCAAAGCTAGAGCTTATGCTGAGTTCTCCAGCCAGCAAAGGAGCGTCTGAGACTTCA
 GCAGCCCTAGCCAGAAGGTTAGTCATCAGTGCCCTTGGGGTTCGAAGTAAAGCAGAGCAAAACCGAACGA
 GAAGCAGAGCTCAAGAACTGCAAGAAGCCAGAGAGAGAAAGCGGTTGGAAGCCAAGCAACGGGAAGAC
 ATCTGGGAAGGCAGAGACCAGTCTACAGTTGA

Restriction Sites: SgfI-MluI

ACCN: NM_001256619

Insert Size: 2379 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001256619.1](#)

RefSeq Size: 3483 bp

RefSeq ORF: 2379 bp

Locus ID: 27291

UniProt ID: [Q7Z5L2](#)

Cytogenetics: 10q24.2

MW: 87.9 kDa