

Product datasheet for **SC330471**

GIDRP88 (R3HCC1L) (NM_001256620) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GIDRP88 (R3HCC1L) (NM_001256620) Human Untagged Clone
Tag:	Tag Free
Symbol:	GIDRP88
Synonyms:	C10orf28; GIDRP86; GIDRP88; PSORT
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC330471 representing NM_001256620.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGCAGCAAGAATCAGAGAGATGCAGAGTTAGAGCCAGAAGGCCTGACATGGCACTTTATGTACCTAAA
GCTCGTAGGGGTGCAGTACTCCTTAAGACAGGTGATGAAGAAGAAAGCTGTGGTTACCTAACTGTG
GTGAAAGAAAAGCAAAAAGAAAGTTCTCTCTCCCAAAAAGAAAGCTTTAAAGACAAACCGGAGGCTCGA
AGACTAAATATCAATCCTGATAGAAAGGAGCATAATTGTAGAGAAGAAAAGAAATCTTCAACAAAATTA
AGAATGGACACATGCCTTCAAAAACAAATAGAGTTTGTCTAAGAGAGGAACCACTGAATCCAAAGAA
GTATTATCCCAAGGACAACAGCAGGGAGCCCCAAATGCTGGGGTTATAACTAATGCACCTTTGCAGAGA
CATTTTAAACCAAGAAGGTGGAGTGTTTGAAGTTGAACTACGGATGTGACAGGACATGAGAGGATA
CTTCTTTACAGGCCTGTTTAAAGATCAGCAGAGGCTCAAGTTCCAAGCAAAACATTCCAAAATGTGGAA
TTCTGTGACTTCAGTAGGCATGAACCTGATGGGAAGCATTGGAAGACAAAGATTGGAAGGCAGAAAT
GAACTGATACCAAGTTTTGGAGATACTATAGAGTTTCTAGAGTTTTAGTTCTGTCATGAAACCT
GAGAATATGATTGACCAATAAACTAAGCTCTGATTCTGAAATTGTACAACAAAGCATGCAACATCA
GATGGAATATTGAATCCAGCAGCGGAGGCATCACCCTACTTCTGTTCTGGAAGTCCAGATGGTGTC
TTTGATCAAACTTGCCTAGATTTTGAAGTTGAGAGTGTAGGTGGTATAGCAATAGTACAGGTTTCATC
TTAGATCAAAAAGATACAGATTCCATTCTGCAACTATGGGTACATCTCTCTGTCAGAGAGCACAAAT
GACACTGTTAGTCCAGTAATGATTAGAGAATGTGAGAAGAATGACAGCACTGCTGATGAGTTACATGTA
AAGCACGAACCTCCTGATACAGCTGTCCTTGCTCATGAAACACATAGAGATAGTGGATTTAAGAATGTA
GGTGACATTACCAATAAAGCATGTATGATGGACACTACAGGTATGTCCTGTAGTGATCATGTAACGTG
GATAGCCCTTATGTAGTTGCAGTTAGAATAGCTGATGAGACCTCTATTAATACACGAAGTTTCTCAAAG
TTTGATAGGAATGAGTGCAGATGCAACCCCTCTTCATGTAGCTAGAAGTGGGAATGACACTGAAGATTTC
AGCAACCCCTTCTGCTTGCTCAGATATTTATGGTGAGAGTATTTTCATCTCATTTTACAGAGTCAACAGGA
AAGTTGATAGAGAGCTTGTGAGATTGTGCTTCTCCTTACCTATAAAAAGATTGCTGGTAGTAATTAT
AACACTTTTTTGGACTCTGAACTCAGTATGTTAAATGGGACAAAAGTTCTTTAGACAGTGCCGTGGGC
ATTGACCTGGGTAGTACTGGTGATACAACAGAAGCATTGCACGAAGTAAAGAACTGCCGAAGATTCAAA
ACAGAAGAGCAAGATGACTCAGGGAGTATAGAATTTGGTGATCTTTTCTGATAGGGAATCATCATCT
ATGGAAACATCCATCGAACCAAAAGCAACTGAACTTCTCACACAGAGGGAATTACTGCCATTGAGGAG
AGCTGGGAGTCTATGTTTAAAGATGATGGTGACTGCCTGGATCCACGTCTTCTACAAGAGTTATCAGGG
AATACCAAGAGCAGAGAGAGCATCCAGGAACCTAGATCTGATTACTACAATCATGAAGTTCTGATATT
GACCTCAGTGATTGTGAATTCACATGTCATTGAAATTTATGACTTTCCCAAGAAATTCGTACTGAA
GACCTTCTACGGGTTTTCTGCAGTTATCAAAAGAAAGGATTTGATATTAATGGGTGGATGATACACAT
GCCCTAGGAGTATTCTCCAGTCCAATTACAGCTCGTGATGCGTTGGGTATTAACACACCATGGTGAAG
ATTCGTCCTTGTACAGGCCACAAGAGCAGCCAAAGGCAAGCTAGAGCTTATGCTGAGTTCCTCCAG
CCAGCAAGGAGCGTCCTGAGACTTCAGCAGCCCTAGCCAGAAGGTTAGTCATCAGTGCCCTTGGGGTT
CGAAGTAAGCAGAGCAAAACCGAACGAGAAGCAGAGCTCAAGAACTGCAAGAAGCCAGAGAGAGAAAG
CGGTTGGAAGCCAAGCAACGGGAAGACATCTGGAAGGCAGAGACCAGTCTACAGTTGA
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001256620

Insert Size: 2337 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001256620.1](#)

RefSeq Size: 3390 bp

RefSeq ORF: 2337 bp

Locus ID: 27291

UniProt ID: [Q7Z5L2](#)

Cytogenetics: 10q24.2

MW: 86.3 kDa