

Product datasheet for **SC330472**

GIDRP88 (R3HCC1L) (NM_001256621) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GIDRP88 (R3HCC1L) (NM_001256621) Human Untagged Clone
Tag:	Tag Free
Symbol:	GIDRP88
Synonyms:	C10orf28; GIDRP86; GIDRP88; PSORT
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC330472 representing NM_001256621. Blue=Insert sequence Red=Cloning site Green=Tag(s)

ATGTTATCAGGGAATACCAAGAGCAGAGAGCATCCAGGAACCTAGATCTGATTACTACAATCATGAA
GTTCTGATATTGACCTCAGTGATTGTGAATCCACATGTCATTGAAATTTATGACTTTCCCAAGAA
TTTCGTAAGACCTTCTACGGGTTTCTGCAGTTATCAAAGAAAGGATTTGATTTAAATGGGTG
GATGATACACATGCCCTAGGAGTATTCTCCAGTCCAATTACAGCTCGTGATGCGTTGGGTATTAACAC
ACCATGGTGAAGATTCGTCCTTGTACAGGCCACAAGAGCAGCAAGGCCAAAGCTAGAGCTTATGCT
GAGTTCTCCAGCCAGCAAAGGAGCGTCTGAGACTTCAGCAGCCCTAGCCAGAAGGTTAGTCATCAGT
GCCCTTGGGGTTGGAAGTAAGCAGAGCAAAACCGAACGAGAAGCAGAGCTCAAGAACTGCAAGAAGCC
AGAGAGAGAAAGCGGTTGGAAGCCAAGCAACGGGAAGACATCTGGGAAGGCAGAGACCAGTCTACAGTT
TGA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001256621
Insert Size:	555 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001256621.1

RefSeq Size: 1587 bp

RefSeq ORF: 555 bp

Locus ID: 27291

UniProt ID: Q7Z5L2

Cytogenetics: 10q24.2

MW: 21.2 kDa