

Product datasheet for **SC330484**

STOML1 (NM_001256677) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: STOML1 (NM_001256677) Human Untagged Clone
Tag: Tag Free
Symbol: STOML1
Synonyms: hUNC-24; SLP-1; STORP
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC330484 representing NM_001256677.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

ATGCCAGATGTACCCAGAGCTGGCCCTCTGCCTCTGTCATGGCCTCATCAGTTTCTGGGGTTCTTG
CTGCTGTTGGTCACCTTCCCATTTCTGGCTGGTTTGCCTGAAGATTGTGCCACCTACGAGCGGATG
ATTGTGTTCCGCCTGGCCGGATCCGCACCCCCAGGGACCTGGCATGGTTCTGCTCTTGGCCTTCATT
GACTCCTTTCAGAGGGTGGATCTGAGGACACGAGCCTTCAACGTCCTCCCTGCAAGCTGGCCTCTAAG
GACGGGGCTGTGCTGTCCGTGGGAGCCGATGTCCAGTTTCGCATCTGGGACCCGGTGTGTCGGTGATG
ACTGTGAAAGACCTGAACACAGCCACACGCATGACAGCCAGAAGCCATGACCAAGGCCCTGCTCAAG
AGGCCGCTGCGGGAGATCCAGATGGAGAAGCTCAAGATCAGCGACCAGCTTCTGCTGGAGATCAACGAT
GTGACCAGGGCCTGGGGGCTGGAGGTAGACCGCTGGAGCTGGCAGTGGAGGCCGTGCTCCAGCCGCC
CAGGACAGCCAGCTGGGCCAACCTGGACAGCACCTCCAGCAGCTGGCCCTGCACTTCTCTGGGAGGA
AGCATGAACTCAATGGCAGGAGGTGCCCGTCCCCGGGGCCAGACACCGTGGAGATGGTGAAGTT
GAGCCACCTGCCCTCAAGTTGGTGCCAGGTCCAGTCCGAAGCAGCCTCTGGCGGAGGGGCTACTGACT
GCTCTACAGCCCTTCTGTCTGAGGCCCTGGTCCAGCAAGTCGGGGCCTGCTACCAAGTTCAATGTCGTC
CTGCCCAGCGGCACCCAAAGCGCCTACTTCTGGACCTCACTACAGGACGAGGAAGAGTGGGACACGGG
GTGCTGATGGCATCCCTGATGTGGTGGTGGAGATGGCCGAGGCAGACCTGCGGGCCCTGCTATGCAGA
GAGCTGCGGCCCTGGGGCCTACATGAGTGGACGGCTGAAGGTGAAGGGCGACCTGGCTATGCCATG
AAGCTGGAGGCTGTCTCAGGGCCTGAAGTAG
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001256677

Insert Size: 1068 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001256677.1](#)

RefSeq Size: 2169 bp

RefSeq ORF: 1068 bp

Locus ID: 9399

UniProt ID: [Q9UBI4](#)

Cytogenetics: 15q24.1

Protein Families: Transmembrane

MW: 38.5 kDa

Gene Summary: May play a role in cholesterol transfer to late endosomes (PubMed:19696025). May play a role in modulating membrane acid-sensing ion channels. Can specifically inhibit proton-gated current of ASIC1 isoform 1. Can increase inactivation speed of ASIC3. May be involved in regulation of proton sensing in dorsal root ganglions (By similarity). May play a role in protecting FBXW7 isoform 3 from degradation (PubMed:23082202).[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (7) contains an alternate exon in the 5' region, initiates translation at an alternate start codon and uses an alternate in-frame splice site in the coding region, compared to variant 1. the encoded isoform (7) is shorter and has a distinct N-terminus, compared to isoform 1.