

Product datasheet for **SC330537**

USP17L13 (NM_001256855) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: USP17L13 (NM_001256855) Human Untagged Clone
Tag: Tag Free
Symbol: USP17L13
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC330537 representing NM_001256855.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

ATGGAGGAGGACTCACTCTACTTGGGTGGTGAGTGGCAGTTCAACCACCTTTCAAACCTCACATCTTCT
CGGCTCGATGCAGCTTTTGCTGAAATCCAGCGGACTTCTCTCCCTGAGAAGTCACCACTCTCATGTGAG
ACCGTGTCGACCTCTGTGATGATTGGTTCCTGAGGCAAGACAGCTTGCTCCCAGGGAGAAGCTTCCT
CTGAGTAGCAGGAGACCTGCTGCGGTGGGGCTGGGCTCCAGAATATGGGAAATACCTGCTACGTGAAC
GCTTCCTTGCAGTGCCTGACATACACACCGCCCTTGCCAACATACATGCTGTCCCGGGAGCACTCAA
ACGTGTCACTCGTACAAGGGCTGCATGCTGTACTATGCAAGCTCACATCACACGGGCCCTCCACAAT
CTGGCCACGTCATCCAGCCCTCACAGGCATTGGCTGCTGGCTTCCATAGAGGCAAGCAGGAAGATGCC
CATGAATTTCTCATGTTCACTGTGGATGCCATGAAAAAGGCATGCCTTCCCGGGCACAAGCAGGTAGAT
CATCCCTCTAAGGACACCACCCTCATCCACCAATATTTGGAGGCTACTGGAGATCTCAAATCAAGTGT
CTCCACTGCCACGGCATTTCAGACACTTTTGACCCTTACCTGGACATCGCCCTGGATATCCAGGCAGCT
CAGAGTGTCAGCAAGCTTTGGAACAGTTGGTGAAGCCCGAAGAAGCTCAATGGAGAGAATGCCTATCAT
TGTGGTGTGTCTCCAGAGGGCGCCGGCCTCCAAGACGTTAACTTTACACACCTCTGCCAAGGTCTC
ATCCTTGATTGAAGAGATTCTCGATGTGACAGGCAACAAGATTGCCAAGAATGTGCAATATCCTGAG
TGCTTGACATGCAGCCATACATGTCTCAGCAGAACACAGGACCTTTGTCTATGCTCTATGCTGTG
CTGGTCCACGCTGGGTGGAGTTGTCAACGACATTACTTCTCTATGTCAAAGCTCAAGAAGGCCAA
TGGTATAAAATGGATGATGCCGAGGTCAACGCCCTAGCATCACTTCTGTCTGAGTCAACAGGCCTAC
GTCTCTTTTACATCCAGAAGAGTGAATGGGAAAGACACAGTGAGAGTGTGTCAAGAGGCAGGGAACCA
AGAGCCCTTGGCGCAGAAGACACAGACAGGCGAGCAACGCAAGGAGAGCTCAAGAGAGACCACCCCTGC
CTCCAGGCCCCGAGTTGGACGAGCACTTGGTGGAAAGAGCCACTCAGGAAAGCACCTTAGACCGCTGG
AAATTCCTTCAAGAGCAAAACAAAACGAAGCCTGAGTTCAACGTGAGAAAGTCGAAGGTACCTGCCT
CCCAGCTACTTGTGATTCAATCAAATACAAGTGTGGGATGAAGAACCATCATCTGAACAGCAA
AGCTCCCTGCTAAACCTCTCTTCGTCGACCCGACACATCAGGAGTCCATGAACACTGGCACACTCGT
TCCCTGCGAGGGAGGCCAGGAGATCCAAGGGAAGAACAAACACAGCAAGAGGGCTCTGCTTGTGTGC
CAGTGA

Restriction Sites: SgfI-MluI
ACCN: NM_001256855
Insert Size: 1593 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001256855.1</u>
RefSeq Size:	1593 bp
RefSeq ORF:	1593 bp
Locus ID:	100287238
UniProt ID:	<u>C9JLJ4</u>
Cytogenetics:	4p16.1
MW:	59.7 kDa
Gene Summary:	Deubiquitinating enzyme that removes conjugated ubiquitin from specific proteins to regulate different cellular processes that may include cell proliferation, progression through the cell cycle, apoptosis, cell migration, and the cellular response to viral infection. [UniProtKB/Swiss-Prot Function]