

## Product datasheet for **SC330560**

### CCDC51 (NM\_001256965) Human Untagged Clone

#### Product data:


**Product Type:** Expression Plasmids  
**Product Name:** CCDC51 (NM\_001256965) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** CCDC51  
**Synonyms:** MITOK  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC330560 representing NM\_001256965.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGGTGGCTCGAGGGCTTGTCAGAGGGCTCGGGAGGACTTGAAGTTCACCAGGCCAAGCTGAAGGAG
GTGAGGGACCGCTTGACCGTGTCTCCAGGGAGGACAGTCAGTACTTGAAGTGGCTACTCTCGAGCAC
AGGATGCTGCAGGAGGAGAAGAGGCTTCGCACAGCCTATCTGCGTGCAGAAGACTCTGAGCGAGAGAAG
TTCTCCCTCTTCTGCAGCTGTGCGGAAAGTCATGAGAAGGAGCGCACAGGGCTGAGAGGACCAAG
AACTGGTCCCTCATTGGCTCAGTCTGGGGCCCTGATTGGTGTGGCTGGCTCCACCTATGTGAACCGT
GTGCGACTACAGGAGCTGAAGGCTTTACTCCTGGAGGCGCAGAAGGGCCTGTGAGTCTCAAGAGGCC
ATTCGAGAACAGGCGTCTAGTACTCCCGCCAGCAGAGGGACCTCCACAATCTCATGGTGGACTTGAGG
GGCCTGGTACATGCTGCTGGGCCAGGGCAGGACTCTGGGTCACAGGCAGGTAGTCCCCGACCAGAGAC
AGAGATGTAGATGTCCTTTCAGCTGCCTTAAAAGAGCAGCTTAGTCATTCCAGGCAAGTCCATTCATGT
CTAGAAGGCTTACGAGAGCAGCTTGATGGCCTAGAAAAGACTTGTAGCCAAATGGCTGGGGTGGTTCAG
CTTGTAAGTCTGCAGCACACCCAGGCCTGGTGAACAGCAGACGGGGCTATGCCAGCTTCTTGCTG
GAGCAGGGGAGCATGATCTTGGCACTGTCAGACACGGAGCAGAGACTAGAAGCCCAAGTCAACAGGAAC
ACCATCTATAGCACCCCTGGTCACCTGTGTGACATTTGTGGCCACACTGCCTGTGCTCTACATGCTATTC
AAAGCCAGCTAA
```

**Restriction Sites:** Sgfl-Mlul



[View online »](#)

**Plasmid Map:**


**ACCN:** NM\_001256965

**Insert Size:** 909 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq:	<u>NM_001256965.1</u>
RefSeq Size:	1611 bp
RefSeq ORF:	909 bp
Locus ID:	79714
UniProt ID:	<u>Q96ER9</u>
Cytogenetics:	3p21.31
Protein Families:	Transmembrane
MW:	33.7 kDa
Gene Summary:	<p>Mitochondrial potassium channel located in the mitochondrial inner membrane (PubMed:31435016). Together with ABCB8/MITOSUR, forms a protein complex localized in the mitochondria that mediates ATP-dependent potassium currents across the inner membrane (that is, mitoK(ATP) channel) (PubMed:31435016). May contribute to the homeostatic control of cellular metabolism under stress conditions by regulating the mitochondrial matrix volume (PubMed:31435016).[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (3) uses two alternate splice sites and initiates translation at a downstream, in-frame start codon, compared to variant 1. Variants 3, 4, 5, 6 and 7 encode the same isoform (2), which has a shorter N-terminus compared to isoform 1.</p>