

Product datasheet for **SC330574**

DDX6 (NM_001257191) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: DDX6 (NM_001257191) Human Untagged Clone
Tag: Tag Free
Symbol: DDX6
Synonyms: HLR2; IDDILF; P54; RCK; Rck/p54
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC330574 representing NM_001257191.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

ATGAGCACGCCAGAACAGAGAACCCTGTTATAATGGGTCTGTCCAGTCAAAATGGTCAGCTGAGAGGC
CCTGTGAAACCCACTGGTGGCCCTGGAGGAGGGGGCACACAGACACAGCAACAGATGAACCACTGAAA
AACACCAACACAATCAATAATGGCACTCAGCAGCAAGCACAGAGTATGACCACCACTATTAACCTGGT
GATGACTGGAAAAAGACTTTAAACTCCCTCAAAGGATCTAAGAATCAAACTTCGGATGTGACCTCC
ACAAAAGGAAATGAGTTGAAGATTACTGTTTGAACGGGAGTTACTGATGGGAATTTTGAATGGGC
TGGGAAAAGCCATCTCCTATTCAGGAGGAGAGCATTCCCATTGCTTTATCTGGTAGGGATATCTTAGCT
AGAGCAAAAAATGGAACAGGCAAGAGCGGTGCCTACCTCATTCCCTTACTTGAACGGCTAGACCTGAAG
AAGGACAATATACAAGCAATGGTGATTGTTCCCACTAGAGAAGTTGCTCTACAGGTCAGTCAAATTTGC
ATCCAGGTGAGCAAAACACATGGGAGGGGCCAAAGTGATGGCAACCACAGGAGGAACCAATTTACGAGAT
GACATAATGAGGCTTGATGATACAGTGCACGTGGTGATTGCTACCCCTGGGAGAAATCCTGGATCTTATT
AAGAAAGGAGTAGCAAAGTTGATCATGTCCAGATGATAGTATTGGATGAGGCAGATAAGTTGCTGTCA
CAGGATTTGTGCAGATAATGGAGGATATTATTCTCACGCTACCTAAAAACAGGCAGATTTTACTATAT
TCCGCTACTTTCCCTCTTAGTGTACAGAAGTTCATGAATCCCATTTCAGAAAACCCTATGAGATTAAC
CTGATGGAGGAACTAAGTCTGAAGGGAGTAACCCAGTACTACGCATATGTAAGTACGCGCCAAAAAGTA
CACTGCCTCAACACACTTTTCTCCAGGCTTCAGATAAACACAGTGCATTTTTCTGTAAGTCTCTCAG
CGAGTTGAATTGCTAGCCAAGAAGATTTCTCAACTGGGTTATTCTTGCTTCTATATTCATGCTAAAATG
AGGCAGGAACATCGAAATCGTGTATTTTCATGATTTCCGAAATGGCTTATGCCGCAATCTTGTGTTGCACT
GATCTGTTTACCCGAGGTATTGATATACAAGCTGTGAATGTGGTAATAAACTTTGATTTCCCAAAGCTG
GCAGAGACCTATCTCCATCGTATTGGAAGATCAGGTCGCTTTGGTCATCTTGGCTTAGCCATCAACTTG
ATCACATATGATGATCGTTCAACCTGAAAAGTATTGAGGAGCAGCTGGGAACAGAAATTAACCTATT
CCGAGCAACATTGATAAGAGCCTGTATGTGCCAGAATACCACAGCGAGCCTGTAGAAGATGAGAAACCT
TAA
  
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001257191

Insert Size: 1452 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001257191.1](#)

RefSeq Size: 6173 bp

RefSeq ORF: 1452 bp

Locus ID: 1656

UniProt ID: [P26196](#)

Cytogenetics: 11q23.3

Protein Pathways: RNA degradation

MW: 54.4 kDa

Gene Summary: This gene encodes a member of the DEAD box protein family. The protein is an RNA helicase found in P-bodies and stress granules, and functions in translation suppression and mRNA degradation. It is required for microRNA-induced gene silencing. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Mar 2012]
Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1 and 2 encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.