

## Product datasheet for **SC330633**

### HARS1 (NM\_001258040) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** HARS1 (NM\_001258040) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** HARS1  
**Synonyms:** CMT2W; HARS; HRS; USH3B  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC330633 representing NM\_001258040.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```


ATGGCAGAGCGTGCGGCGCTGGAGGAGCTGGTGAACCTCAGGGAGAGCGGTGCGAGGCCCTCAAGCAG
CAGAAGGCCAGCGCCGAGCTGATCGAGGAGGAGGTGGCGAACTCCTGAACTGAAGGCACAGCTGGGT
CCTGATGAAAGCAAACAGAAATTTGTGCTCAAACCCCAAGGAAACTGATGGGAAAGTATGGGGAA
GACTCCAAGCTTATCTATGACCTGAAGGACCAGGGCGGGAGCTCCTGTCCCTTCGCTATGACCTCACT
GTTCTTTTGTGCGTATTTGGCAATGAATAAACTGACCAACATTAACGCTACCACATAGCAAAGGTA
TATCGGCGGGATAAACCAGCCATGACCCGTGGCCGATACCGGGAATTCTACCAGTGTGATTTTGACATT
GCTGGGAACTTTGATCCCATGATCCCTGATGCAGAGTGCCTGAAGATCATGTGCGAGATCCTGAGTTCA
CTTCAGATAGGCGACTTCTGGTCAAGGTAACGATCGACGCATTCTAGATGGGATGTTTGTCTATCTGT
GGTGTTCCTGACAGCAAGTTCCTGACCATCTGCTCCTCAGTAGACAAGCTGGACAAGGTGTCCTGGGAA
GAGGTGAAGAATGAGATGGTGGGAGAGAAGGGCCTTGACCTGAGGTGGCTGACCGCATTGGGGACTAT
GTCCAGCAACATGGTGGGTATCCCTGGTGAACAGCTGCTCCAGGATCCTAAACTATCCAAAACAAG
CAGGCCTTGGAGGGCCTGGGAGACCTGAAGTTGCTCTTTGAGTACCTGACCCTATTTGGCATTGATGAC
AAAATCTCCTTTGACCTGAGCCTTGCTCGAGGGCTGGATTACTACACTGGGGTATCTATGAGGCAGTG
CTGCTACAGACCCAGCCAGGCAGGGGAAGGCCCTGGGTGTGGCAGTGTGGCTGCTGGAGGACGC
TATGATGGGCTAGTGGCATGTTTCGACCCAAAGGGCGCAAGGTGCCATGTGTGGGCTCAGCATTGGG
GTGGAGCGGATTTTCTCCATCGTGAACAGAGACTAGAGGCTTTGGAGGAGAAGTACGGACCACGGAG
ACACAGGTGCTTGTGGCATCTGCACAGAAGAAGCTGCTAGAGGAAAGACTAAAGCTTGTCTCAGAACTG
TGGGATGCTGGGATCAAGGCTGAGCTGTGTACAAGAAGAACCCTAAAGCTACTGAACCAGTTACAGTAC
TGTGAGGAGGCAGGCATCCCACTGGTGGCTATCATCGGCGAGCAGGAACTCAAGGATGGGGTATCAAG
CTCCGTTCAGTGACGAGCAGGGAAGAGGTGGATGTCCGAAGAGAAGACCTTGTGGAGGAAATCAAAGG
AGAACAGGCCAGCCCTCTGCATCTGTA
  
```

**Restriction Sites:** Sgfl-Mlul



[View online »](#)

## Plasmid Map:



ACCN: NM\_001258040

Insert Size: 1410 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001258040.1](#)

RefSeq Size: 2202 bp

RefSeq ORF: 1410 bp

Locus ID: 3035

UniProt ID: [P12081](#)

**Cytogenetics:** 5q31.3

**Protein Pathways:** Aminoacyl-tRNA biosynthesis

**MW:** 52.7 kDa

**Gene Summary:** Aminoacyl-tRNA synthetases are a class of enzymes that charge tRNAs with their cognate amino acids. The protein encoded by this gene is a cytoplasmic enzyme which belongs to the class II family of aminoacyl-tRNA synthetases. The enzyme is responsible for the synthesis of histidyl-transfer RNA, which is essential for the incorporation of histidine into proteins. The gene is located in a head-to-head orientation with HARSL on chromosome five, where the homologous genes share a bidirectional promoter. The gene product is a frequent target of autoantibodies in the human autoimmune disease polymyositis/dermatomyositis. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2012]

Transcript Variant: This variant (2) lacks an alternate in-frame exon compared to variant 1. The resulting isoform (2) has the same N- and C-termini but is shorter compared to isoform 1.