

Product datasheet for **SC330649**

WDR89 (NM_001258272) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	WDR89
Synonyms:	C14orf150; MSTP050
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC330649 representing NM_001258272. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGAAGATTGAGGAACAATTTGCTAATCTGCACATTGTTAAATGTTCTTAGGAACCAAGAGCCC
ACTTACCTTCTTGGTATAGACACATCAAAGACTGTCCAAGCAGGAAAGGAAACTTGGTTGCTGTTTAA
TGTTCTAATGGATCAATCAGAATATATGATAAGAAAGGTTAAATGTACTACGAGAATTTAGTGGATAT
CCTGGACTTCTTAATGGAGTCAGATTTGCAAATTCCTGTGACAGTGATATTCAGCATGTACTGATGGC
ACTGTGAAATGCTGGGATGCTCGAGTAGCCAGAGAAAAACCTGTTCAAGGTTACCCCTTCC
AATATTTTTATCAGTTTTGATATTAATTGTAATGATCATATTATTTGTGCTGGTACAGAAAAAGTTGAT
GATGATGCATTGTTGGTGTGTTGGGATGCAAGGATGAATTCAGAAATTTATCTACAATAAGACTCA
CTTGGTGCAATTCAGAGACATAGTGATGATGTCCTCAAGTACGTTTCCATCCCAGCAATCCCAAC
ATGGTAGTCTCAGGTTTCATCTGATGGCCTGGTAAATGATTTGATATTAATATTGATAATGAGGAGGAT
GCACTGGTTACAACCTGTAACCTCAATTTTCATCAGTAAGCTGTATTGGTGGTCTGGGAAAGGTTATAAA
CAGATTTACTGCATGACACATGATGAAGGATTTTATTGGTGGGATCTTAATCATCTGGACACTGATGAA
CCAGTTACACGTTTGAACATCCAGGATGTCAGAGAAGTAGTTAACATGAAAGAAGATGCTTTGGACTAT
TTGATTGGTGGCCTATATCATGAAAAGACAGACATTGCATGTTATTGGAGGAACAAACAAAGGAAGG
ATTCATTTGATGAAGTGCAGCATGTCAGGACTGACCCATGTGACTAGCCTTCAGGGAGGGCATGCTGCT
ACAGTCCGTTCTTTCTGTTGGAATGTGCAAGATGATTCTTTGTTGACTGGAGGAGAAGATGCACAGTTG
TTACTTTGGAACCTGGAGCTATAGAGAAGACCTTTACAAAGAAAGAGAGTATGAAAAATGCATCCTCT
GTGCACCAACGAGTACGAGTTCATAGTAATGATTCTTATAAAAGAAGGAAAAAGCAGTGA

```

Restriction Sites:	SgfI-MluI
ACCN:	NM_001258272
Insert Size:	1164 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001258272.1
RefSeq Size:	3680 bp
RefSeq ORF:	1164 bp
Locus ID:	112840
UniProt ID:	Q96FK6
Cytogenetics:	14q23.2
MW:	43.2 kDa