

## Product datasheet for **SC330666**

### SLC44A3 (NM\_001258340) Human Untagged Clone

#### Product data:

|               |   |
|---------------|---|
| Product Type: | Expression Plasmids                         |
| Product Name: | SLC44A3 (NM_001258340) Human Untagged Clone |
| Tag:          | Tag Free                                    |
| Symbol:       | SLC44A3                                     |
| Synonyms:     | CTL3  |
| Vector:       | pCMV6-Entry (PS100001)                      |




[View online »](#)

Fully Sequenced ORF: >SC330666 representing NM\_001258340.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGCACTGCCTGGGCGCCGAGTACCTGGTTTCTGCAGAAGGAGCCCTAGGCAAAGGGAGTGGCGACCC
CAGATTTATAGGAAATGCACAGATACGGCATGGTTATTCTGTCTTTCTCTTTGGACTGGTTTGGTG
TTTATCATGGGCTACTCGGTGGTGGCTGGAGCCGCGGAAGACTCCTCTTTGGCTATGACAGCTTTGGC
AACATGTGTGGCAAGAAGAAGAACTCCCCCGTGAAGGGGCCCTCTTTCAGGGCAGGACATGACCCTAAAA
AAACACGTGTTCTTTATGAATTCCTGCAACCTGGAAGTCAAAGGTACGCAGCTCAACCGCATGGCCCTC
TGTGTATCCAACCTGCCCTGAAGAGCAGCTTGACTCCCTGGAAGAGGTCCAGTTCTTTGCAAAACACAGT
GGGTCCTTCTGTGTGTTTATAGTTTGAATTCCTTCAACTATACCCACAGTCCAAAAGCAGACTCACTG
TGTCCCAGGCTACCAGTTCCTCCAAGCAAGTCAATTCCTTATTTAACCGATGTGTCCTCAAACACCT
GAGTGCTACTCCCTATTTGCATCTGTTTTGATAAATGATGTTGACACCCTCCACCGAATTCTAAGTGGA
ATCATGTCGGGAAGAGATACAATCCTTGGCCTGTGTATCCTCGCATTAGCCTTGCTTTGGCCATGATG
TTTACCTTCAGATTCATACCACCCTTCTGGTTCACATTTTCATTTTATTGGTTATTTGGGATTGTTG
TTTGTCTGCGGTGTTTTATGGTGGCTGTATTACTATACCAACGACCTCAGCATAGAATTGGACACA
GAAAGGGAAAATATGAAGTGCCTGCTGGGGTTTGTATCGTATCCACAGGCATCAGGCAGTGTGCTC
GTCTTGATTTTTGTTCTCAGAAAAGAAATAAAATTGACAGTTGAGCTTTTCCAAATCACAAATAAAGCC
ATCAGCAGTGTCCCTTCTGCTGTTCCAGCCACTGTGGACATTTGCCATCCTCATTTTCTTCTGGGTC
CTCTGGGTGGCTGTGCTGCTGAGCCTGGGAAGTGCAGGAGCTGCCCAGGTATGGAAGGCGGCCAAGTG
GAATATAAGCCCTTTCGGGCATTCGGTACATGTGGTTCGTACCATTTAATTGGCCTCATCTGGACTAGT
GAATTCATCCTTGCCTGCCAGCAAATGACTATAGCTGGGCAGTGGTACTTGTATTTCACAGAAGT
AAAAATGATCCTCCTGATCATCCCATCCTTTCGTCTCTCCTCCTCTCTTCTTCTACCATCAAGGAACC
GTTGTGAAAGGGTCATTTTTAATCTCTGTGGTGGGATTCCGAGAATCATTGTCATGACATGCAAAAC
GCACTGAAAGAACAGCATGGTGCAATTGTCCAGGTACCTGTTCCGATGCTGCTACTGCTGTTTCTGGTGT
CTTGACAAATACCTGCTCCATCTCAACCAGAATGCATATACTACAACCTGCTATTAATGGGACAGATTTT
TGTACATCAGCAAAGATGCATTCAAAATCTTGTCCAAGAAGTCAAGTCACTTTACATCTATTAAGTGC
TTTGGAGACTTCATAATTTTTCTAGGAAAGGTGTAGTGGTGTGTTTCACTGTTTTTGGAGGACTCATG
GCTTTTAACTACAATCGGGCATTCCAGGTGTGGCAGTCCCTCTGTTATTGGTAGCTTTTTTGCCTAC
TTAGTAGCCCATAGTTTTTATCTGTGTTTGAAGTGTGCTGGATGCACCTTTCTGTGTTTTGCTGTT
GATCTGGAACAAATGATGGATCGTCAGAAAAGCCCTACTTTATGGATCAAGAATTTCTGAGTTTCGTA
AAAAGGAGCAACAAATTAACAATGCAAGGGCACAGCAGGACAAGCACTCATTAAAGGAATGAGGAGGGA
ACAGAACTCCAGGCCATTGTGAGTAG
```

Restriction Sites: Sgfl-Mlul

**Plasmid Map:**


**ACCN:** NM\_001258340

**Insert Size:** 1959 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001258340.1](#)

**RefSeq Size:** 2394 bp

**RefSeq ORF:** 1959 bp

**Locus ID:** 126969

**UniProt ID:** [Q8N4M1](#)

**Cytogenetics:** 1p21.3

**Protein Families:** Transmembrane

**MW:** 73.7 kDa