

Product datasheet for **SC330667**

SLC44A3 (NM_001258341) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: SLC44A3 (NM_001258341) Human Untagged Clone
Tag: Tag Free
Symbol: SLC44A3
Synonyms: CTL3
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC330667 representing NM_001258341.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```


ATGCACTGCCTGGGCGCCGAGTACCTGGTTTCTGCAGAAGGAGCCCTAGGCAAAGGGAGTGGCGACCC
CAGATTTATAGGAAATGCACAGATACGGCATGGTTATTCTGTTCTTTCTCTTTGGACTGGTTTGGTG
TTTATCATGGGCTACTCGGTGGTGGCTGGAGCCGCGGAAGACTCCTCTTTGGCTATGACAGCTTTGGC
AACATGTGTGGCAAGAAGAACTCCCCGTGGAAGGGGCCCTCTTCAGGGCAGGACATGACCCTAAAA
AAACACGTGTTCTTTATGAATTCCTGCAACCTGGAAGTCAAAGGGTCTTCTCTGTGTGTTATAGTTG
AATTCCTCAACTATACCCACAGTCCAAAAGCAGACTCACTGTGTCCAGGCTACCAGTTCCTCCAAGC
AAGTCATTTCCCTTATTTAACCGATGTGTCCCTCAAACACTGAGTGCTACTCCCTATTTGCATCTGTT
TTGATAAATGATGTTGACACCCTCCACCGAATTCTAAGTGAATCATGTGCGGAAGAGATACAATCCTT
GGCCTGTGTATCCTCGCATTAGCCTTGTCTTTGGCCATGATGTTACCTCAGATTCATCACCACCCTT
CTGGTTCACATTTTCATTTTCATTGGTTATTTGGGATTGTTGTTGTCTGCGGTGTTTTATGGTGGCTG
TATTATGACTATACCAACGACCTCAGCATAGAATTGGACACAGAAAGGGAAAATAGAAGTGCCTGCTG
GGGTTTGCTATCGTATCCACAGGCATCACGGCAGTGCTGCTGCTTGTATTTTGTCTCAGAAAGAGA
ATAAAATTGACAGTTGAGCTTTTCCAAATCACAATAAAGCCATCAGCAGTGTCCCTTCTGCTGTTT
CAGCCACTGTGGACATTTGCCATCCTCATTTTCTCTGGGTCTCTGGGTGGCTGTGCTGAGCCTG
GGAAGTGCAGGAGCTGCCAGGTTATGGAAGGCGCCAAGTGAATATAAGCCCTTTTCGGGCATTTCGG
TACATGTGGTTCGTACCATTTAATTGGCCTCATCTGGACTAGTGAATTCATCCTTGCCTGCCAGAAATG
ACTATAGCTGGGGCAGTGGTTACTTGTATTTCAACAGAAGTAAAAATGATCCTCCTGATCATCCCATC
CTTTCGTCTCTCCATTCTCTTCTTCTACCATCAAGGAACCGTTGTGAAAGGGTCATTTTAACTCT
GTGGTGAGGATCCGAGAATCATTGTACATGACATGCAAAACGCACTGAAAGAACAGCATGGTGCATTG
TCCAGGTACCTGTCCGATGCTGCTACTGCTGTTTCTGGTGTCTTGACAAATACCTGCTCCATCTCAAC
CAGAATGCATATACTACAAGTCTATTAATGGGACAGATTTCTGTACATCAGCAAAAGATGCATTCAA
ATCTTGTCCAAGAAGTCAAGTCACTTTACATCTATTAAGTCTTTGGAGACTTCATAATTTTCTAGGA
AAGGTGTTAGTGGTGTGTTTCACTGTTTTTGGAGGACTCATGGCTTTTAACTACAATCGGCATTCCAG
GTGTGGGCAGTCCCTCTGTTATTGGTAGCTTTTTTGCCTACTTAGTAGCCCATAGTTTTTATCTGTG
TTTGAAACTGTGCTGGATGCACTTTTCTGTGTTTGTGTTGATCTGGAAACAAATGATGGATCGTCA
GAAAAGCCCTACTTTATGGATCAAGAAATTTCTGAGTTTCGTAAAAAGGAGCAACAAATTAACAATGCA
AGGGCACAGCAGGACAAGCACTCATTAAAGGAATGAGGAGGGAACAGAAGTCCAGGCCATTGTGAGATAG
  
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001258341

Insert Size: 1863 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001258341.1](#)

RefSeq Size: 2298 bp

RefSeq ORF: 1863 bp

Locus ID: 126969

UniProt ID: [Q8N4M1](#)

Cytogenetics: 1p21.3

Protein Families: Transmembrane

MW: 70.1 kDa