

## Product datasheet for **SC330668**

### SLC44A3 (NM\_001258342) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** SLC44A3 (NM\_001258342) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** SLC44A3  
**Synonyms:** CTL3  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC330668 representing NM\_001258342.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```


ATGCACTGCCTGGGCGCCGAGTACCTGGTGTATCATGGGCTACTCGGTGGTGGCTGGAGCCGCGGGA
AGACTCCTCTTTGGCTATGACAGCTTTGGCAACATGTGTGGCAAGAAGAACTCCCCGTGGAAGGGGCC
CCTCTTTAGGGCAGGACATGACCCTAAAAAACACGTGTTCTTTATGAATTCCTGCAACCTGGAAGTC
AAAGGTACGCAGCTCAACCGCATGGCCCTCTGTGTATCCAACCTGCCCTGAAGAGCAGCTTGACTCCCTG
GAAGAGGTCCAGTCTTTGCAAACACCAGTGGTCCCTCTGTGTGTTTATAGTTTGAATTCCTTCAAC
TATACCCACAGTCCAAAAGCAGACTCAGTGTGCCAGGCTACCAGTTCCTCAAGCAAGTCATTTCCC
TTATTTAACCGATGTGTCCCTCAAACACCTGAGTGCTACTCCCTATTTGCATCTGTTTTGATAAATGAT
GTTGACACCCTCCACCGAATTCTAAGTGAATCATGTCCGGAAGAGATACAATCCTTGGCCTGTGTATC
CTCGCATTAGCCTTGTCTTTGGCCATGATGTTTACCTTCAGATTCATCACCACCCTTCTGGTTACATT
TTCATTTTCATTGGTTATTTTGGGATTGTTGTTGTCTGCGGTGTTTTATGGTGGCTGATTATGACTAT
ACCAACGACCTCAGCATAGAATTGGACACAGAAAGGAAAATATGAAGTGCCTGCTGGGTTTGCTATC
GTATCCACAGGCATCACGGCAGTGTCTGCTGCTGATTTTTGTTCTCAGAAAGAGAATAAAATTGACA
GTTGAGCTTTTCAAATCACAAATAAAGCCATCAGCAGTGTCTCCCTTCTGCTGTTCCAGCCACTGTGG
ACATTTGCCATCCTCATTTTCTTGGTCTCTGGGTGGCTGTGCTGCTGAGCCTGGAACCTGCAGGA
GCTGCCAGGTTATGGAAGGCGCCAAGTGAATATAAGCCCTTTCGGGCATTCCGTACATGTGGTCCG
TACCATTTAATTGGCCTCATCTGGACTAGTGAATTCATCCTTGCCTGCCAGCAAATGACTATAGCTGGG
GCAGTGGTTACTTGTTATTTCAACAGAAGTAAAAATGATCCTCCTGATCATCCCACCTTTCGTCTCTC
TCCATTCTTCTTCTACCATCAAGGAACCGTTGTGAAAGGGTCATTTTTAATCTCTGTGGTGAGGATT
CCGAGAATCATTGTCATGTACATGCAAAACGCACTGAAAGAACAGCAGCATGGTGCATTGTCAGGTAC
CTGTTCCGATGCTGCTACTGCTGTTTCTGGTGTCTTGACAAATACCTGCTCCATCTCAACCAGAATGCA
TATACTACAACCTGCTATTAATGGGACAGATTTCTGTACATCAGCAAAAGATGCATTCAAAATCTTGCC
AAGAACTCAAGTCACTTTACATCTATTAAGTCTTTGGAGACTTCAATAATTTTTCTAGGAAAGGTGTTA
GTGGTGTGTTTCACTGTTTTTGGAGGACTCATGGCTTTTAACTACAATCGGGCATTCCAGGTGTGGCA
GTCCCTCTGTTATTGGTAGCTTTTTTGGCTACTTAGTAGCCCATAGTTTTTATCTGTGTTTGAACCT
GTGCTGGATGCACTTTTCTGTGTTTGTGCTGTTGATCTGGAAACAAATGATGGATCGTCAGAAAAGCCC
TACTTTATGGATCAAGAATTTCTGAGTTTCGTAAGGAGCAACAAATTAACAATGCAAGGGCACAG
CAGGACAAGCACTCATTAAGGAATGAGGAGGGAACAGAAGTCCAGGCCATTGTGAGATAG
  
```



[View online »](#)

Restriction Sites: Sgfl-MluI

Plasmid Map:



ACCN: NM\_001258342

Insert Size: 1854 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001258342.1](#)

RefSeq Size: 2289 bp

RefSeq ORF: 1854 bp

Locus ID: 126969

UniProt ID: [Q8N4M1](#)

Cytogenetics: 1p21.3

Protein Families: Transmembrane

MW: 69.4 kDa

